

# Vasc BAND<sup>™</sup> hemostat



The clear difference in hemostasis


# Vasc BAND™ hemostat

## Inflatable compression

**Detachable strap**  
secures syringe during  
patient transport

**Extended hook & loop retention straps**  
provide quick application and  
enhance patient comfort

**Clear plastic band**  
provides access  
site visualization



**Locking plunger syringe**  
enables precise 2ml deflation increments


**Single compression balloon**  
delivers optimal hemostasis  
pressure for every patient

**Secure twist connector**  
provides easy and reliable connection  
between syringe and inflation valve

# for access site hemostasis

**Five band sizes**  
accommodate any wrist size:  
Short - 21cm  
Regular - 24cm  
Long - 27cm  
Extra Long - 29cm  
XXL - 37cm

The Vasc Band delivers more direct radial compression than the standard inflatable radial compression device using recommended clinical protocols, with minimal compression of the ulnar artery.\*




**Test protocol:** The Vasc Band and the standard inflatable radial compression devices were tested on 10 human subjects. For each test, the compression band was secured to the wrist and inflated with recommended nominal and maximum inflation values for each device. Per the product IFU's, nominal air inflation is 15ml for the Vasc Band and 13ml for the standard inflatable radial compression device, while maximum air inflation of 18ml is the same for both devices. Pressure applied by each device to the radial and ulnar arteries was recorded using a Tekscan ELF Force and Load Measurement System and FlexiForce B201 flexible sensors. Data shown in the chart reflects maximum inflation value of 18ml for both devices.

\*Data reference: SR1700 (VSI Vasc Band, Competitive Product Performance Study)

# Vasc BAND hemostat

The Vasc Band hemostat is a compression device to assist hemostasis of arterial, venous and hemodialysis percutaneous access sites.

Model	Band Size	Box Quantity
3521	Short - 21cm	10
3524	Regular - 24cm	10
3527	Long - 27cm	10
3529	Extra Long - 29cm	10
3537	XXL - 37cm	10



Please see the Instructions for Use for a complete listing of the indications, contraindications, warnings and precautions.

CAUTION: Federal Law (U.S.A.) restricts this device to sale by or on the order of a physician.

Vasc Band is a trademark of Vascular Solutions, Inc.

©2016 Vascular Solutions, Inc. All rights reserved. ML2752 Rev. H 11/16

This literature is intended for distribution only in the U.S. and select international markets.



Manufactured for  
**Vascular Solutions, Inc.**  
6464 Sycamore Court North  
Minneapolis, Minnesota 55369 USA

Customer Service:  
United States: 888.240.6001  
customerservice@vasc.com  
www.vasc.com