



LMA® Airway Portfolio

Inspiring confidence right from the start







300m

LMA<sup>®</sup> Airways have been used in 300 million patients since 1988

Every 3 seconds, an LMA<sup>®</sup> Airway is used somewhere in the world\*



LMA<sup>®</sup> is the most extensively researched laryngeal mask brand in the world

# Seals, silicone and safety features: The latest innovations in LMA® Airway design

## It starts with a high-quality seal

It's one of the most critical aspects of laryngeal mask performance – a reliable oropharyngeal seal to help protect a patient's airway during general anesthesia. That's why the new products in the LMA® Airway portfolio from Teleflex incorporate the latest innovations in cuff design – features designed to help you establish, verify, and maintain a seal you can count on.

## Silicone:

### A material that gently conforms to the anatomy

The original laryngeal mask airways were made of silicone – a material valued for its softness and flexibility.<sup>1,2</sup>

Teleflex recognizes the benefits of silicone and has introduced silicone products into the single-use LMA® Airway portfolio.

Potential benefits of silicone over polyvinyl chloride (PVC) include:

- Greater elasticity to conform to the anatomy<sup>1-3</sup>
- Higher oropharyngeal seal pressures<sup>2,4</sup>
- A lower risk of post-operative sore throat<sup>5,6</sup>
- Not made with phthalates

**100%**  
**Silicone Cuff**





**Cuff Pilot™ Technology:  
Because intracuff pressure counts**


An optimal intracuff pressure of 40-60 cm H<sub>2</sub>O is important to achieve a reliable oropharyngeal seal with a laryngeal mask. Over-inflation of a cuff can have a negative impact on oropharyngeal seal pressures and may put patients at risk of a wide range of airway morbidities.<sup>7-10</sup> Studies have shown laryngeal mask intracuff pressures exceed 60 cm H<sub>2</sub>O in a vast majority of patients.<sup>11-15</sup>

Cuff Pilot™ Technology was developed to support clinicians in avoiding the known risks of cuff hyperinflation, which include sore throat, dysphagia, an increased risk of aspiration due to leakage around the cuff, and hypoglossal, lingual, or recurrent laryngeal nerve palsies.<sup>7-10</sup> Cuff Pilot™ Technology is an integrated, single-use cuff pressure indicator that constantly monitors cuff pressure, detecting changes resulting from fluctuations in temperature, nitrous oxide levels, and movements within the airway.

**Simple to use:** a simple glance is all that's needed to confirm optimal intracuff pressure right from the start – and throughout the procedure

**Supports a reliable seal:** designed to help clinicians maintain the optimal intracuff pressure of 40-60 cm H<sub>2</sub>O required to achieve a reliable oropharyngeal seal

**Strategy to reduce the risk of airway morbidity:** throughout a procedure clinicians can monitor intracuff pressure to ensure it does not exceed 60 cm H<sub>2</sub>O – a strategy shown to reduce airway morbidity by up to 70%<sup>9-11</sup>




# LMA<sup>®</sup> Gastro<sup>™</sup> Airway

## Proactive airway management for endoscopy procedures

Hypoxemia due to respiratory depression or airway obstruction is a known risk associated with endoscopic procedures, with studies showing that hypoxemia can occur in 11–50% of cases.<sup>16-18</sup> Currently, many endoscopy procedures are undertaken without an airway management device in place; therefore, when hypoxemia occurs, failure to rescue the airway can lead to potential complications and undesirable outcomes.<sup>19</sup>

The LMA<sup>®</sup> Gastro<sup>™</sup> Airway with Cuff Pilot<sup>™</sup> Technology from Teleflex is the first laryngeal mask specifically designed to enable clinicians to proactively manage their patients' airways while facilitating direct endoscopic access via the integrated endoscope channel. Once general anesthesia is established, the laryngeal mask can be placed and the endoscope passed through the large bore endoscope channel, providing direct access to the esophagus and upper gastrointestinal tract.

Cuff Pilot<sup>™</sup> Technology



Endoscope channel

Adjustable holder and strap


Integrated bite block

Silicone airway tube and cuff

# LMA® Unique™ (Silicone Cuff) Airway

## A first-generation airway with the added benefits of integrated cuff pressure monitoring


The LMA® Unique™ (Silicone Cuff) Airway with Cuff Pilot™ Technology from Teleflex combines all the features of the LMA® Unique™ Airway with the added benefits of integrated cuff pressure monitoring and a silicone cuff.<sup>1,2,14</sup> Cuff Pilot™ Technology was developed to support clinicians in avoiding the known risks of cuff hyperinflation, which include sore throat, dysphagia, an increased risk of aspiration due to leakage around the cuff, and hypoglossal, lingual or recurrent laryngeal nerve palsies.<sup>7-10</sup> Silicone cuffs are designed to better conform to patient anatomy, compared with PVC cuffs, which may result in improved airway sealing characteristics.<sup>2</sup> The LMA® Unique™ (Silicone Cuff) Airway is an ideal choice for routine anesthetic procedures, for difficult airway situations, or for airway management during cardiopulmonary resuscitation.



# LMA® Unique EVO™ Airway

## Versatile airway for routine use, with intubation capabilities for difficult airway situations

The LMA® Unique™ EVO Airway with Cuff Pilot™ Technology from Teleflex is designed to provide the versatility needed to adapt to changing clinical situations. As a first-generation laryngeal mask, the LMA® Unique™ EVO Airway is suitable for routine use, while providing fiberscope-assisted intubation capabilities to give clinicians additional flexibility to manage the airway during unforeseen difficult airway situations. A fixed-curve airway tube facilitates rapid insertion, while the silicone cuff is soft and flexible, and conforms to the anatomy to create an effective oropharyngeal seal.<sup>1,2</sup>



# LMA® Protector™ Airway

## Expanding applications with functional separation of the respiratory and GI tracts


The LMA® Protector™ Airway with Cuff Pilot™ Technology from Teleflex is the most advanced laryngeal mask in the LMA® Airway portfolio. It features an elongated silicone cuff designed to form an optimal dual seal to functionally separate the respiratory tract from the gastrointestinal tract. The LMA® Protector™ Airway is the first single-use, second-generation laryngeal mask to feature a dual gastric drainage channel and pharyngeal chamber intended to channel gastric contents away from the airway to help reduce the risk of aspiration.<sup>20</sup> The LMA® Protector™ Airway also offers fiberscope-assisted direct intubation capabilities for additional flexibility to manage the airway during unexpected difficult airway situations, and is ideal for routine use or more challenging procedures.



# Innovation, performance, reliability

You'll find it all right here

Thirty years after the introduction of the first laryngeal mask airway, Teleflex continues to lead the way. With millions of uses around the world, supported by thousands of clinical studies, the LMA® Airway portfolio is proven to address a broad range of challenges across multiple clinical situations. Moving forward, we will continue our relentless pursuit of identifying unmet clinical needs just as we have right from the start.



# Ordering information

## LMA<sup>®</sup> Gastro<sup>™</sup> Airway with Cuff Pilot<sup>™</sup> Technology\*

ITEM NUMBER	MASK SIZE	PATIENT WEIGHT (KG)	MAXIMUM INTRACUFF PRESSURE (CM H <sub>2</sub> O*)	MAXIMUM ENDOSCOPE SIZE (OD) (MM)
1E5030	3	30-50	60	14
1E5040	4	50-70	60	14
1E5050	5	70-100	60	14

## LMA<sup>®</sup> Unique<sup>™</sup> (Silicone Cuff) Airway with Cuff Pilot<sup>™</sup> Technology\*

ITEM NUMBER	MASK SIZE	PATIENT WEIGHT (KG)	MAXIMUM INTRACUFF PRESSURE (CM H <sub>2</sub> O)
105200-000010	1	Up to 5	60
105200-000015	1.5	5-10	60
105200-000020	2	10-20	60
105200-000025	2.5	20-30	60
105200-000030	3	30-50	60
105200-000040	4	50-70	60
105200-000050	5	70-100	60
105200-000060	6	More than 100	60

## LMA<sup>®</sup> Unique EVO<sup>™</sup> Airway with Cuff Pilot<sup>™</sup> Technology\*

ITEM NUMBER	MASK SIZE	PATIENT WEIGHT (KG)	MAXIMUM INTRACUFF PRESSURE (CM H <sub>2</sub> O)	MAXIMUM ETT ID (MM)
1D2030	3	30-50	60	7
1D2040	4	50-70	60	7.5
1D2050	5	70-100	60	8

## LMA<sup>®</sup> Protector<sup>™</sup> Airway with Cuff Pilot<sup>™</sup> Technology\*

ITEM NUMBER	MASK SIZE	PATIENT WEIGHT (KG)	MAXIMUM INTRACUFF PRESSURE (CM H <sub>2</sub> O)	MAXIMUM ETT ID (MM)	LARGEST SIZE OG TUBE (FR)
192030	3	30-50	60	6.5	16
192040	4	50-70	60	7.5	18
192050	5	70-100	60	7.5	18

ETT=Endotracheal Tube OG=Orogastric Tube

\* Cuff Pilot<sup>™</sup> Technology recommendation for cuff pressure: green zone pressure range = 40-60 cm H<sub>2</sub>O for all sizes.



[teleflex.com/lma](https://teleflex.com/lma)

## References:

1. Eschertzhuber S, Brimacombe J,† Hohlrieder M, Keller C.‡ *Anaesthesia*. 2009;64(1):79-83.
2. Jagannathan N, Sohn LE, Sawardekar A, Gordon J, Langen KE, Anderson K. *Anaesthesia*. 2012;67(6):632-639.‡
3. Anand LK, Goel N, Singh M, Kapoor D. *Acta anaesthesiologica Taiwanica*. 2016;54(2):44-50.
4. Maitra S, Khanna P, Baidya DK. *Eur J Anaesthesiol*. 2014;31(5):266-273.
5. Wong JG, Heaney M, Chambers NA, Erb TO, von Ungern-Sternberg BS. *Paediatr Anaesth*. 2009;19(5):464-469.
6. Bell SF, Morris NG, Rao A, Wilkes AR, Goodwin N. *Anaesthesia*. 2012;67(12):1337-1342.
7. van Zundert T, Brimacombe J.‡ *Anaesthesiol Intensive Ther*. 2012;44(2):63-70.‡
8. Burgard G, Mollhoff T, Prien T. *J Clin Anesth*. 1996;8(3):198-201.
9. Chantzara G, Stroumpoulis K, Alexandrou N, Kokkinos L, Iacovidou N, Xanthos T. *Minerva Anesthesiol*. 2014;80(5):547-555.
10. Seet E, Yousaf F, Gupta S, Subramanyam R, Wong DT, Chung F. *Anesthesiology*. 2010;112(3):652-657.‡
11. Vasanth Karthik R, Ranganathan P, Kulkarni AP, Sharma KS. *J Anesth*. 2014;28(5):662-667.
12. Schloss B, Rice J, Tobias JD. *Int J Pediatr Otorhinolaryngol*. 2012;76(2):284-286.
13. von Ungern-Sternberg BS, Erb TO, Chambers NA, Heaney M. *Paediatr Anaesth*. 2009;19(9):837-843.
14. Martin DP, Bhalla T, Thung A, Tobias JD. *Int J Pediatr Otorhinolaryngol*. 2013;77(1):76-78.‡
15. Rokamp KZ, Secher NH, Moller AM, Nielsen HB. *BMC Anesthesiol*. 2010;10:20.
16. Amornyotin S. Cardiorespiratory complications during moderate and deep sedation for gastrointestinal endoscopic procedures. In: Amornyotin S, ed: InTech; 2013: <https://www.intechopen.com/books/endoscopy-of-gi-tract/cardiorespiratory-complications-during-moderate-and-deep-sedation-for-gastrointestinal-endoscopic-pr>. Accessed 10 July 2017.
17. Cote GA, Hovis RM, Ansstas MA, et al. *Clin Gastroenterol Hepatol*. 2010;8(2):137-142.
18. de Paulo GA, Martins FP, Macedo EP, Goncalves ME, Mourao CA, Ferrari AP. *Endosc Int Open*. 2015;3(1):E7-E13.
19. Qadeer MA, Rocio Lopez A, Dumot JA, Vargo JJ. *Dig Dis Sci*. 2009;54(5):1035-1040.
20. Sng BL, Ithnin FB, Mathur D, Lew E, Han NR, Sia AT. *BMC Anesthesiol*. 2017;17(1):26.‡

† Formerly a paid consultant of Teleflex or its affiliates.

‡ Research sponsored in part by Teleflex Incorporated.

§ Research sponsored by Teleflex Incorporated.

Teleflex is a global provider of medical technologies designed to improve the health and quality of people's lives. We apply purpose driven innovation – a relentless pursuit of identifying unmet clinical needs – to benefit patients and healthcare providers. Our portfolio is diverse, with solutions in the fields of vascular and interventional access, surgical, anesthesia, cardiac care, urology, emergency medicine and respiratory care. Teleflex employees worldwide are united in the understanding that what we do every day makes a difference. For more information, please visit [teleflex.com](http://teleflex.com).

Teleflex is the home of Arrow®, Deknatel®, Hudson RCI®, LMA®, Pilling®, Rüschi®, and Weck® – trusted brands united by a common sense of purpose.

## Corporate Office

Phone +1 610 225 6800, 550 E. Swedesford Road, Suite 400, Wayne, PA 19087, USA

For more information, please visit [teleflex.com](http://teleflex.com).

**Caution:** Federal law (USA) restricts this device to sale by or on the order of a physician.

Teleflex, the Teleflex logo, Arrow, Cuff Pilot, Deknatel, Hudson RCI, LMA, LMA Classic, LMA Gastro, LMA Fastrach, LMA Flexible, LMA ProSeal, LMA Protector, LMA Supreme, LMA Unique, Piling, Rüschi, Unique EVO, and Weck are trademarks or registered trademarks of Teleflex Incorporated or its affiliates in the U.S. and/or other countries. Information in this material is not a substitute for the product Instructions for Use. Not all products may be available in all countries. Please contact your local representative. All data current at time of printing (07/2018). Subject to technical changes without further notice.

© 2018 Teleflex Incorporated. All rights reserved. MC-003980