




RÜSCH EZ-BLOCKER™ KIT

Instructions for Use	en
Gebrauchsanweisung	de
Mode d'emploi	fr
Instrucciones de uso	es
Istruzioni per l'uso	it
Instruções de utilização	pt
Bruksanvisning	sv
Gebruiksaanwijzing	nl
Kullanım Talimatları	tr
Instrukcja obsługi	pl
Инструкция по применению	ru
取扱説明書	ja
使用说明	zh

RÜSCH EZ-BLOCKER™ KIT

en

Instructions for Use

The EZ-Blocker kit contains the following sterile products:

- 1 EZ-Blocker
- 1 EZ-Multiport™ Adapter with scope lid
- 1 Oxygen Adapter
- 2 Dust caps

Caution:

Federal (USA) law restricts this device to sale by or on the order of a physician.

Description of the Device:

The EZ-Blocker is an endobronchial blocking device intended for use in airway management.

The EZ-Blocker is an endobronchial blocker with two distal extensions, both with an inflatable cuff and a central lumen. The shaft, the two distal extensions and the two distal tips are made of radiopaque material. The two distal extensions are colored differently for identification purposes. In addition "blue stripes" are used as symbols to identify the blue distal extension for CPAP purposes (Fig. 2). These color and symbolic identifications (Fig. 2) are repeated on each corresponding balloon assembly for cuff inflation. Printed depth markers indicate the distance to the distal tip of the EZ-Blocker. This device is intended for use through a standard single lumen endotracheal tube to provide lung isolation. Fiberoptic or video bronchoscopy should be used for placement and proper positioning. The EZ-Multiport Adapter is designed for connection to ventilation devices, and allows introduction of a fiberoptic or video bronchoscope and the EZ-Blocker, or a suction catheter.

Indications for Use:

The EZ-Blocker is indicated for use to intubate the patients bronchi in order to differentially isolate the left or right lung for procedures which require one-lung ventilation, lung separation. Patient population: Patients requiring one lung isolation.

Environment of use: Hospitals – OR Operating Room and ICU Intensive Care Unit.

Contraindications:

- Internal endotracheal tube diameter too small to facilitate both the EZ-Blocker and a fiberoptic or video bronchoscope.
- Airway diameter insufficient to allow introduction and passage of the EZ-Blocker.
- When equipment for fiberoptic or video bronchoscope is not available.
- Use of LASER or electro-surgical electrode in the immediate area of the EZ-Blocker.
- Known allergy to polyurethane.

Warnings:

- When both devices are situated in the tube at the same time this may result in hypoxemia and high ventilation pressures.
- The surgeon should be aware of the location of the cuff in order to prevent damage, especially while working near the hilum.
- In the case of a sudden drop in saturation, or rise in ventilation pressure, immediate fiberoptic or video bronchoscope control is mandatory.
- The EZ-Blocker™ should be advanced with care and never be forced down the tube.
- The EZ-Blocker should not be advanced if any resistance is encountered.
- If ventilation worsens acutely during endobronchial blockade, the cuff should be deflated immediately and the cuff position should be inspected.
- Adequate ventilation may be impaired due to the presence of a fiberoptic or video bronchoscope and the EZ-Blocker in the endobronchial tube, at the same time. This may result in higher peak- and expiratory pressures.

- This device is intended for single use only. Do not reuse, re-sterilise or autoclave. Reuse of single use devices creates a potential risk of patient or user infections. Contamination of the device may lead to injury, illness or death of the patient. Cleaning, disinfection and sterilization may compromise essential material and design characteristics leading to device failure. The manufacturer will not be responsible for any direct, incidental or consequential damages resulting from resterilisation or reuse.
- No lidocaine or lubricant containing lidocaine should be used on product balloons.

Precautions:

- It is recommended to use pulse oxymetry when using the device.
- During procedures near the hilum caution is recommended. The position of the cuff should be verified by a surgeon to prevent inadvertent cuff damage.
- During the use of this device for one-lung ventilation, an increased FIO₂ for the ventilated lung is often needed in order to maintain adequate arterial oxygen saturation. In some circumstances however, an increased FIO₂ may not be sufficient to maintain adequate arterial oxygen saturation and other techniques, e.g. CPAP through the central lumen of the EZ-Blocker, may be necessary.
- Care should be taken to ensure the cuff remains fully inflated during longer procedures.
- The patient should be paralyzed during one-lung ventilation.
- The initial placement of the EZ-Blocker should be with the patient in a supine position.
- The position of the inflated cuff should be reconfirmed after positioning and after repositioning the patient for the procedure, especially after turning the patients head.
- The position of the endotracheal tube should be secured with tape before and during the use of the EZ-Blocker.
- During use of the EZ-Blocker the capnograph waveform may become distorted.
- The EZ-Blocker can be introduced through a standard endotracheal tube. The largest size possible fitting the patient is recommended for optimal performance.
- A tube size 7 (inner diameter in mm) is the smallest recommended endotracheal tube size when using the EZ-Blocker.
- The outer diameter of the bronchoscope should be 3,2 to 4,2 mm.

These recommendations are only designed to serve as a general guideline. They are not intended to supersede institutional protocols or professional clinical judgments concerning patient care.

Instructions for Use:

- Take the EZ-Blocker™ and remove EZ-Blocker from cardboard

ATTENTION:

Prior to use remove the protection tube as follows:

- Hold bifurcation in one hand and remove protection tube by pulling on warning label, check EZ-Blocker for damage.
- Check both cuffs for leakage and asymmetric inflation by injecting air into the cuffs with a clean syringe. Place the caps to seal off both proximal endings (luer lock fitting) of the lumina in the EZ-Blocker (Fig. 2).
- Warning: to prevent airway damage, the balloon and cuff should never be overinflated.
- Care must be taken to avoid damage of the product and its accessories.
- After intubation of the patient with a proper sized endotracheal tube, the EZ-Multiport™ Adapter (Fig. 1) should be placed on the endotracheal tube (port D, Fig. 4) and the patient should be ventilated (port A, Fig. 4).
- Deflate both cuffs of the EZ-Blocker completely.
- Lubricate the distal part of the EZ-Blocker and the bronchoscope according to standard practice procedures.
- No lidocaine or lubricant containing lidocaine should be used on product balloons.
- During introduction of the blocking device, **the endotracheal tube should be positioned 4 cm above the carina to ensure proper functioning of the EZ-Blocker** (Fig. 3).
- Introduce the EZ-Blocker through port (B), after removal of the closing plug of port (B), of the EZ-Multiport Adapter and advance into the endotracheal tube (Fig. 5).
- Introduce the fiberoptic or video bronchoscope through the other port (C) of the EZ-Multiport Adapter after removing the plug of the closing cap in order to visualize the airway and the EZ-Blocker (Fig. 6).
- Advance the EZ-Blocker under direct visual control until both extensions are just outside the endotracheal tube (Fig. 7).

- By checking and adjusting the EZ-Blocker position under direct visual guidance, advance the device until both distal extensions have been introduced in both main stem bronchi (Fig. 8).
- Following correct placement under direct visual guidance, inflate the appropriate distal extension cuff by inflating the corresponding balloon; distal cuff on blue extension corresponds with blue striped proximal balloon (Fig. 2). Isolate the target lung for lung isolation (Fig. 9). Auscultation will confirm adequate ventilation of the contra lateral lung.
- Repeated warning: to prevent airway damage, the balloon and cuff should never be overinflated.
- After reaching the appropriate position for the procedure inspect for proper position, signs of cuff herniation and air leakage.
- After repositioning the patient or the patient's head, check the cuff position with the fiberoptic or video bronchoscope again.
- Once the EZ-Blocker™ is in the correct position in the airway of the patient, the cap, mounted on the shaft of the EZ-Blocker (Fig. 5) can be tightened on the port of the EZ-Multiport™ Adapter until air tight.
- It is recommended to make a final check on the position of the inflated cuff, after which the fiberoptic or video bronchoscope is removed and surgery can start.
- In case of repositioning the patient, the closing cap needs to be unscrewed to the extent that the EZ-Blocker can move freely through the closing cap.
- In case of ventilation problems or failure, deflate the cuff and reposition under direct visual guidance.
- Once one lung ventilation is not required any longer, deflate the cuffs completely and remove the blocker from the airway.

Packaging:

The EZ-Blocker™ set is sterilized in ethylene oxide and supplied in peel open package. The entire set is intended for single use only! Do not use if package has been previously opened or damaged. Do not re-sterilize. Do not re-use and do not use after the expiratory date indicated on the package.

Storage Conditions:

The EZ-Blocker set shall be stored in a room suited for storage and at a temperature not exceeding +25°C.

Transport Conditions:

During transport of the EZ-Blocker a temperature in the tolerance range between +40°C and -10°C must be maintained. However, the products may temporarily (not longer than 72 hours) be exposed to temperatures up to +50°C and -29°C.

Rüsch, EZ-Blocker and EZ-Multiport are trademarks or registered trademarks of Teleflex Incorporated or its affiliates. © 2012 Teleflex Incorporated.
All information is in accordance with the state of the art when going to press.
Subject to technical alterations. Printed on chlorine free paper.

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Teleflex Medical Europe Limited
Contract Medical International GmbH