



LMA®
Cuff Pressure Monitoring




The role of laryngeal masks with inflatable cuffs

- Laryngeal masks, including the family of LMA[®] Airways, are routinely used to facilitate oxygenation and ventilation during general anesthesia
 - Laryngeal masks are also increasingly used for emergency airway management in the pre-hospital setting
- Optimal use of a laryngeal mask relies on a number of factors, including selecting a device that is the correct size for a given patient, ensuring that the device is inserted to a correct depth, and achieving an adequate seal between the device and the airway anatomy¹
 - The majority of laryngeal masks include an elliptical cuff that is inflated with air following insertion
 - The inflatable cuff encircles the laryngeal inlet, effectively isolating the distal airways
 - Once inflated, the cuff functions to prevent air from leaking to the atmosphere (Figure 1)

- Optimal inflation of the cuff is essential for patient safety
 - Intracuff pressure must be high enough to seal the airway during both spontaneous and assisted ventilation, but be low enough to avoid reducing/occluding blood flow in the laryngopharyngeal mucosa and/or damaging adjacent nerves

Figure 1. *In situ* positioning of a laryngeal mask with an inflatable cuff [LMA[®] Unique[™] (Silicone Cuff) Airway with Cuff Pilot[™] Technology] depicted



The complex relationship between cuff inflation volume and cuff pressure

- Typically, the inflation system of a laryngeal mask comprises a cuff, an inflation line, a pilot balloon (which provides an indication of the pressure within the cuff), and a check valve (which prevents leakage of air and maintains intracuff pressure)
 - To inflate the cuff, a given volume of air is injected via the inflation line
 - The cuff should be inflated with sufficient air to obtain a low-pressure seal
- Laryngeal mask manufacturers provide recommendations about safe maximum inflation volumes for their devices; however, these are based on the physical properties of the cuff (i.e., the volume to which the cuff can be safely distended without compromising the material)²
 - Typically, the recommended (i.e., maximum) filling volume is employed in clinical practice, even though it is not an indication of what is suitable for the majority of patients²
 - An *in vitro* experiment that used pediatric-sized single-use and reusable laryngeal masks showed that, when starting from a completely deflated cuff, inflation to the maximum recommended volume almost always resulted in an intracuff pressure that was higher than recommended (i.e., >60 cm H₂O)³
 - ♦ Furthermore, when starting from a resting cuff (i.e., with the pilot balloon valve opened to atmospheric pressure), inflation to the maximum recommended volume resulted in a cuff pressure of >120 cm H₂O for all but one of the devices studied³
- An *in vivo* study that used single-use and reusable laryngeal masks in pediatric patients showed similar results; when starting from a completely deflated cuff, the recommended intracuff pressure (i.e., 60 cm H₂O) was exceeded “well below” the recommended maximum inflation volume⁴
 - ♦ Indeed, an intracuff pressure of 60 cm H₂O was achieved with approximately one-half of the recommended maximum inflation volume⁴
- An intracuff pressure of 60 cm H₂O and higher may have clinical consequences, including increased leakage around the cuff^{5,6}
 - Excessive intracuff pressures may also lead to pharyngolaryngeal morbidity, including post-operative sore throat, dysphagia, dysphonia, and/or nerve injury⁷⁻¹⁵



The prevalence of cuff hyperinflation

Table 1. Prevalence of cuff hyperinflation (intracuff pressure ≥ 60 cm H₂O) in children and adults whose airway was managed with various laryngeal masks (comparative data for endotracheal tubes [ETTs] have been included, where available)

STUDY	SETTING	INSERTION/INFLATION METHOD	MEASUREMENT OF CUFF PRESSURE	DEVICE	INTRACUFF PRESSURE	RATE OF CUFF HYPERINFLATION IN %
Infants and children						
von Ungern-Sternberg BS, et al. 2009 ²¹	Elective surgery requiring general anesthesia	Inserted unchanged straight from the sterile packaging without further inflation or deflation of the cuff	Measured using a calibrated hand-held manometer following device insertion	LMA® Classic™ Airway (n=87) LMA® Unique™ Airway (n=89) LMA® Flexible™ Airway [single use] (n=115) LMA Flexible Airway [reusable] (n=80) LMA® ProSeal™ Airway (n=61) PROACT Medical Ltd PRO-Breathe® (n=568)	NR	≥ 60 cm H ₂ O: 21 (67% for size 1 devices)
Schloss B, et al. 2012 ¹⁸	Surgery requiring general anesthesia	Inserted with the cuff partially inflated, as per routine clinical practice, with further inflation as needed to ensure a seal during positive pressure ventilation to a peak inflating pressure of 20–25 cm H ₂ O	Measured using a hand-held manometer within the first 30 minutes	Ambu® A/S laryngeal mask (n=200)	Mean \pm SD: 57 \pm 30	≥ 60 cm H ₂ O: 53
Martin DP, et al. 2013 ²²	Surgery requiring general anesthesia	Inserted with the cuff partially inflated, as per routine clinical practice, with further inflation as needed to ensure a seal during positive pressure ventilation to a peak inflating pressure of 20–25 cm H ₂ O	Measured using a hand-held manometer immediately after device placement	AES Inc laryngeal mask (n=100)	NR	>60 cm H ₂ O: 31
Adults						
Rokamp KZ, et al. 2010 ¹⁶	Elective surgery requiring general anesthesia	Inflated as per the disposition of the head anesthiologist (without the use of a manometer or a pressure release valve)	Measured using a cuff pressure gauge following placement of the airway	Ambu® A/S laryngeal mask (n=82)	Median (range): 95 (10–121)	>60 cm H ₂ O: 68 >120 cm H ₂ O:† 41
				ETT (n=119)	Median (range): 30 (8–100)	>30 cm H ₂ O: 45 >40 cm H ₂ O: 28
Spiro M, et al. 2010 ¹⁹	NR	NR	Measured in the operating room	Fannin Ltd single-use laryngeal mask (n=89)	Median: 120	≥ 60 cm H ₂ O: 85
Patient population not specified						
Sandhu G, et al. 2012 ¹⁷	Surgery requiring general anesthesia	NR	Measured using a hand-held manometer within 30 minutes of device placement	Laryngeal mask (n=34)	NR	≥ 60 cm H ₂ O: 97 60–120 cm H ₂ O: 24 >120 cm H ₂ O:† 74
				ETT (n=27)	Mean: 39	≥ 25 cm H ₂ O: 52
Viernes DC, Joffe AM, et al. 2012 ²⁰	Surgery requiring general anesthesia	NR	Measured using a Compass Lumbar Puncture	Laryngeal masks (n=44)	Median (range): 90 (12–199)	>60 cm H ₂ O: 68 >120 cm H ₂ O: 30
				ETT (n=246)	Median (range): 43 (6–199)	>30 cm H ₂ O: 61 >60 cm H ₂ O: 23

† The upper limit of the pressure gauge/manometer
NR, not reported; SD, standard deviation

- Cuff hyperinflation may be highly prevalent in the setting of general anesthesia
 - Studies have shown that as many as 97% of patients whose airway was managed with a laryngeal mask had a cuff pressure that exceeded the recommended value of no more than 60 cm H₂O (Table 1)^{16–22}
 - Of concern, up to 74% of measurements were more than twice the recommended pressure (i.e., >120 cm H₂O)^{16,17,20}
- Cuff hyperinflation is not solely restricted to the use of laryngeal masks, but is also apparent in patients whose airway is managed with a cuffed endotracheal tube (ETT) (Table 1)^{16,17,20}

The association between intracuff pressure and airway morbidity

Table 2. Frequency of post-operative pharyngolaryngeal complications according to intracuff pressure in patients whose airway was managed with various laryngeal masks

STUDY	SETTING	LARYNGEAL MASK	INSERTION AND INFLATION METHOD	RATE OF AIRWAY-RELATED ADVERSE EVENTS IN %		
Burgard G, et al. 1996 ⁸	Gynecological surgery requiring general anesthesia	LMA [®] Classic™ Airway	Device inserted and inflated with recommended volume of air (25 mL [size 3] or 35 mL [size 4]) Low-pressure group: Intracuff pressure released to minimal airtightness pressure High-pressure group: Intracuff pressure not released		Low-pressure group (n=100)	High-pressure group (n=100)
				Sore throat (minimal, moderate, or severe) in the recovery room and 4, 8, and 24 hours post-surgery	0, 0, 0, and 0	8, 8, 8, and 5
Nott MR, et al. 1998 ¹²	Elective surgery requiring general anesthesia	LMA Classic Airway	Inserted using a standard technique (with a slight lateral approach if resistance encountered) Inflation to move the device into the correct position within the pharynx, up to typical volumes for each size Adjusted group: Intracuff pressure adjusted after 5 minutes until there was a slight leak to positive pressure at 10–12 cm H ₂ O (i.e., “just airtight”) Non-adjusted group: Intracuff pressure left unchanged		Adjusted group (n=412)	Non-adjusted group (n=427)
				Sore throat (mild, moderate, or severe)	7*	16
Seet E, et al. 2010 ¹³	Short-duration elective ambulatory surgery requiring general anesthesia	LMA Classic Airway	Inserted according to the anesthesiologist’s preferred technique and the manufacturer’s instructions Inflated at the discretion of the attending anesthesiologist to achieve an audible seal Pressure-limiting group: Intracuff pressure reduced to 54–60 cm H ₂ O if >60 cm H ₂ O Routine care group: Intracuff pressure left unchanged		Pressure-limiting group (n=97)	Routine care group (n=103)
				Any pharyngolaryngeal complications	13*	46
				Sore throat at 1, 2, and 24 hours post-surgery	7, 2***, and 3**	8, 9, and 14
				Dysphagia at 1, 2, and 24 hours post-surgery	1*, 0*, and 2***	13, 13, and %
Chantzara G, et al. 2014 ⁹	Elective urological surgery requiring general anesthesia	LMA [®] laryngeal mask (distributed by Iamex SA in Greece)	Gradually inflated at the discretion of the anesthesiologist to achieve a seal without audible leak during positive pressure ventilation with a maximum intracuff volume of 30 mL (size 4) and 40 mL (size 5) Intervention group: Intracuff pressure maintained at 60 cm H ₂ O Observation group: Intracuff pressure left unchanged		Intervention group (n=60)	Observation group (n=60)
				Any pharyngolaryngeal adverse effects 24 hours post-surgery	8*	35
				Sore throat at 1, 2, and 24 hours post-surgery	7, 2***, and 3**	8, 9, and 14
				Dysphagia at 1, 2, and 24 hours post-surgery	1*, 0*, and 2***	13, 13, and %
Kang JE, et al. 2014 ¹⁰	Laparoscopic surgery requiring general anesthesia	LMA [®] Supreme™ Airway	Device inserted with the cuff completely deflated Cuff inflated thereafter Low-pressure group: Intracuff pressure limited to 25 cm H ₂ O High-pressure group: Intracuff pressure set to 60 cm H ₂ O		Low-pressure group (n=49)	High-pressure group (n=52)
				Sore throat on days 1 and 2	4 and 6***	12 and 23
				Dysphagia on days 1 and 2	0*** and 0***	8 and 8
				Dysphonia on days 1 and 2	0 and 0	0 and 2
Vasanth Karthik R, et al. 2014 ¹⁵	Short-duration elective surgery requiring general anesthesia	LMA [®] ProSeal™ Airway	Inserted according to the anesthesiologist’s preferred technique Inflated to no more than the maximum recommended volume to achieve a seal without audible leak during positive pressure ventilation with a tidal volume of 8 mL/kg and a peak inspiratory pressure <25 cm H ₂ O Pressure-monitored group: Intracuff pressure reduced to 60 cm H ₂ O if >60 cm H ₂ O Control group: Intracuff pressure left unchanged		Pressure-monitored group (n=60)	Control group (n=59)
				Any pharyngolaryngeal complications	32	42
				Sore throat at 1, 2, and 24 hours post-surgery	5, 22, and 25	9, 37, and 41
				Dysphagia at 1, 2, and 24 hours post-surgery	3, 12, and 17	8, 20, and 22

- There is a large body of evidence to show that excessive intracuff pressures can have a detrimental effect on a patient’s airway
 - Morbidity may manifest as post-operative sore throat, dysphagia, dysphonia, and/or nerve injury⁷⁻¹⁵
- The pharyngeal mucosal perfusion pressure may be exceeded during the use of laryngeal masks²³ and this might lead to cuff pressure-related airway morbidity²⁴⁻²⁶
 - If intracuff pressure exceeds perfusion pressure, the mucosa may become ischemic, leading to tissue damage²⁴
- Other types of pressure-related morbidity include cranial nerve injuries (e.g., the lingual, laryngeal, hypoglossal, and glossopharyngeal nerves)^{7,11,14}
 - Such injuries are thought to be the result of pressure neuropraxia, with hyperinflation of the cuff being a contributing factor^{7,11,14}
- Several randomized controlled trials have been undertaken to study the effect of laryngeal mask cuff pressure on the incidence of post-operative pharyngolaryngeal complaints
 - Data show that a reduction in intracuff pressure results in a decreased rate of airway morbidity (Table 2)^{8-10,12,13,15}
 - In one study,¹³ reducing intracuff pressure to 54–60 cm H₂O led to a 70% reduction in pharyngolaryngeal complications
- A reduction in the rate of airway morbidity when intracuff pressure is reduced to the recommended maximum has also been observed with cuffed ETTs^{27,28}
- Interestingly, rates of sore throat appear to have increased over time, from 13% when the LMA Classic Airway was first described in 1983²⁹ to rates approaching 50% in recent years^{15,30}

Cuff pressure monitoring in clinical practice

- In a clinical setting, intracuff pressure is typically estimated via digital palpation of the pilot balloon
 - Data show that the palpation technique is inaccurate³¹ and tends to result in an underestimation of actual intracuff pressure^{31,32}
- Clinical endpoints (e.g., appropriate positioning of the device and an adequate seal) may also be used to guide cuff inflation
 - This method is associated with cuff hyperinflation in a majority of patients^{6,33} and with increased leakage around the cuff⁶
- In order to avoid cuff hyperinflation, numerous researchers have concluded that intracuff pressure should be routinely monitored/controlled in both adults and children,^{3,4,6,8,13,16,19-21,32-34} typically with a pressure manometer
 - Cuff pressure monitoring is considered particularly important in children
 - ♦ Because pediatric patients have smaller airway diameters than adults, the effect of swelling of the airway is known to be greater^{4,6}
 - ♦ Furthermore, with smaller-sized devices (i.e., those used in pediatric patients), small changes in volume can result in large changes in pressure⁴

Table 3. Changes in cuff pressure in certain surgical settings/scenarios

SURGICAL SETTING/SCENARIO	CHANGE IN CUFF PRESSURE
Nitrous oxide anesthesia	The use of nitrous oxide is known to increase intracuff volume and pressure over time, owing to the more rapid diffusion of nitrous oxide (versus air) across the wall of the cuff ^{35,36†}
Changes in patient position	Significant increases in cuff pressure have been observed following rotation of the head during surgery ³⁷
Changes in atmospheric pressure	Increases in cuff pressure have been observed following increases in altitude/elevation, ³⁸ which is relevant during aeromedical transport

- Because changes in intracuff pressure over time are not uncommon, particularly in certain surgical settings/scenarios (Table 3), monitoring should occur throughout the use of a laryngeal mask^{4,8,16,33}
 - ♦ Additional costs related to the repair, replacement, and maintenance (e.g., calibration) of such devices would also be incurred, as would costs associated with cleaning the devices between each patient^{18,22}
- Despite calls for the use of manometry to monitor intracuff pressure, this is not routinely undertaken in many institutions^{2-4,6,13,18,21,32,33}
 - Medical personnel may consider the use of a manometer time consuming, difficult/cumbersome, and/or inaccurate⁴⁰⁻⁴³
- The lack of uptake of manometers is likely multifactorial
 - Single measurements of intracuff pressure using external manometers do not provide a continuous assessment of cuff pressure throughout the course of surgery³⁰
 - Researchers in the United States note that the cost of a commercially available manometer ranges from approximately \$US100 to \$US400;^{13,18,22} in Europe, the cost is estimated to be approximately €100 per unit¹⁶
 - ♦ Based on these prices, the installation of manometers into every operating room would be costly³⁹ and potentially infeasible^{18,22}

Cuff Pilot™ Technology An integrated pressure indicator

- Effective cuff inflation is about pressure, not volume
- Cuff Pilot™ Technology is an integrated pressure indicator that constantly monitors intracuff pressure and provides at-a-glance manometry of cuff pressure levels
 - Those LMA Airways that include Cuff Pilot Technology can be identified by the inclusion of Cuff Pilot Technology in the product name
- Cuff Pilot Technology replaces the standard pilot balloon and is being introduced on all single-use LMA® Airways with a silicone cuff
 - This is achieved using a color-coded scale in which pressure ranges are indicated using specific color zones (Figure 2)




Figure 2. Integrated Cuff Pilot Technology with a color-coded scale to indicate intracuff pressure



† Diffusion of nitrous oxide is known to differ according to the physical properties of the cuff

Cuff Pilot™ Technology

An integrated pressure indicator

- A number of studies[†] conducted in pediatric and adult patients have evaluated the accuracy, clinical performance, and tolerability of laryngeal masks that include the integrated Cuff Pilot™ Technology^{22,30,44}
 - One study noted that use of Cuff Pilot Technology allowed for an accurate indication of intracuff pressure²²
 - Another study showed that use of Cuff Pilot Technology resulted in a lower incidence of post-operative pharyngolaryngeal complications, compared with a device without an integrated pressure indicator³⁰
- Cuff Pilot Technology can help clinicians avoid many of the potential disadvantages of manometry
 - The use of LMA® Airways with Cuff Pilot Technology will help avoid the need for expenditure related to the purchase, maintenance, repair, cleaning, and storage of manometers, which could result in substantial cost savings
 - The single-use nature of LMA Airways with Cuff Pilot Technology may help reduce the risk of cross-contamination
 - Use of Cuff Pilot Technology helps avoid issues associated with
- the management, inventory, availability, and calibration of manometers within hospital operating rooms
 - Because Cuff Pilot Technology is integrated into all single-use LMA Airways with a silicone cuff, it does not require users to take additional steps to determine cuff pressure
 - Cuff Pilot Technology constantly monitors intracuff cuff pressure, as opposed to a manometer, which provides a spot measurement

References

1. van Zundert AA, Kumar CM, van Zundert TC. Malpositioning of supraglottic airway devices: preventive and corrective strategies. *Br J Anaesth.* 2016;116(5):579-582.
2. Bick E, Bailes I, Patel A, Brain AI. Fewer sore throats and a better seal: why routine manometry for laryngeal mask airways must become the standard of care. *Anaesthesia.* 2014;69(12):1304-1308.
3. Maino P, Dullenkopf A, Keller C, Bernet-Buettiker V, Weiss M. Cuff filling volumes and pressures in pediatric laryngeal mask airways. *Paediatr Anaesth.* 2006;16(1):25-30.
4. Wallace CJ, Chambers NA, Erb TO, von Ungern-Sternberg BS. Pressure volume curves of paediatric laryngeal mask airways. *Anaesthesia.* 2009;64(5):527-531.
5. Hockings L, Heaney M, Chambers NA, Erb TO, von Ungern-Sternberg BS. Reduced air leakage by adjusting the cuff pressure in pediatric laryngeal mask airways during spontaneous ventilation. *Paediatr Anaesth.* 2010;20(4):313-317.
6. Licina A, Chambers NA, Hullett B, Erb TO, von Ungern-Sternberg BS. Lower cuff pressures improve the seal of pediatric laryngeal mask airways. *Paediatr Anaesth.* 2008;18(10):952-956.
7. Bruce IA, Ellis R, Kay NJ. Nerve injury and the laryngeal mask airway. *J Laryngol Otol.* 2004;118(11):899-901.
8. Burgard G, Mollhoff T, Prien T. The effect of laryngeal mask cuff pressure on postoperative sore throat incidence. *J Clin Anesth.* 1996;8(3):198-201.
9. Chantzara G, Stroumpoulis K, Alexandrou N, Kokkinos L, Iacovidou N, Xanthos T. Influence of LMA cuff pressure on the incidence of pharyngolaryngeal adverse effects and evaluation of the use of manometry during different ventilation modes: a randomized clinical trial. *Minerva Anestesiol.* 2014;80(5):547-555.
10. Kang JE, Oh CS, Choi JW, Son IS, Kim SH. Postoperative pharyngolaryngeal adverse events with laryngeal mask airway (LMA Supreme) in laparoscopic surgical procedures with cuff pressure limiting 25 cmH(2)O: prospective, blind, and randomised study. *Sci World J.* 2014;2014:709801.
11. Kawauchi Y, Nakazawa K, Ishibashi S, Kaneko Y, Ishikawa S, Makita K. Unilateral recurrent laryngeal nerve neuropraxia following placement of a ProSeal laryngeal mask airway in a patient with CREST syndrome. *Acta Anaesthesiol Scand.* 2005;49(4):576-578.
12. Nott MR, Noble PD, Parmar M. Reducing the incidence of sore throat with the laryngeal mask airway. *Eur J Anaesthesiol.* 1998;15(2):153-157.
13. Seet E, Yousaf F, Gupta S, Subramanyam R, Wong DT, Chung F. Use of manometry for laryngeal mask airway reduces postoperative pharyngolaryngeal adverse events: a prospective, randomized trial. *Anesthesiology.* 2010;112(3):652-657.
14. Thiruvankatarajan V, Van Wijk RM, Rajbhoj A. Cranial nerve injuries with supraglottic airway devices: a systematic review of published case reports and series. *Anaesthesia.* 2015;70(3):344-359.
15. Vasanth Karthik R, Ranganathan P, Kulkarni AP, Sharma KS. Does cuff pressure monitoring reduce postoperative pharyngolaryngeal adverse events after LMA-ProSeal insertion? A parallel group randomised trial. *J Anesth.* 2014;28(5):662-667.
16. Rokamp KZ, Secher NH, Moller AM, Nielsen HB. Tracheal tube and laryngeal mask cuff pressure during anaesthesia - mandatory monitoring is in need. *BMC Anesthesiol.* 2010;10:20.
17. Sandhu G, Thanawala V, Flack J. Prospective audit to determine endotracheal tube and laryngeal mask airway cuff pressures during general anaesthesia (abstr. 87). *Anaesthesia.* 2012;67(Suppl. s2):67.
18. Schloss B, Rice J, Tobias JD. The laryngeal mask in infants and children: what is the cuff pressure? *Int J Pediatr Otorhinolaryngol.* 2012;76(2):284-286.
19. Spiro M, Gross J, Boomers O. The influence of laryngeal mask airway (LMA) cuff pressure on postoperative sore throat (abstr. 19AP3-1). *Eur J Anaesthesiol.* 2010;27(47):250.
20. Viernes DC, Joffe AM. Routine intracuff monitoring of airway devices should be standard practice (abstr. A578). Abstract presented at the Annual Meeting of the American Society of Anesthesiologists, Washington, DC, USA, 13-17 October, 2012. <http://www.asaabstracts.com/strands/asaabstracts/abstract.htm?sessionid=A5556126EC5F6BCDB5F7D97499F3F268?absnum=4875&index=8&year=2012>. Accessed 15 February, 2016.
21. von Ungern-Sternberg BS, Erb TO, Chambers NA, Heaney M. Laryngeal mask airways--to inflate or to deflate after insertion? *Paediatr Anaesth.* 2009;19(9):837-843.
22. Martin DP, Bhalla T, Thung A, Tobias JD. Clinical evaluation of a novel LMA with a color-coded pressure gauge. *Int J Pediatr Otorhinolaryngol.* 2013;77(1):76-78.
23. Marjot R. Pressure exerted by the laryngeal mask airway cuff upon the pharyngeal mucosa. *Br J Anaesth.* 1993;70(1):25-29.
24. Aydogmus MT, Eksioğlu B, Oba S, et al. Comparison of laryngeal mask airway Supreme and laryngeal mask airway ProSeal for laryngopharyngeal trauma and postoperative morbidity in children. *Braz J Anesthesiol.* 2013;63(6):445-449.
25. Brimacombe J, Holyoake L, Keller C, et al. Pharyngolaryngeal, neck, and jaw discomfort after anesthesia with the face mask and laryngeal mask airway at high and low cuff volumes in males and females. *Anesthesiology.* 2000;93(1):26-31.
26. O'Kelly SW, Heath KJ, Lawes EG. A study of laryngeal mask inflation. Pressures exerted on the pharynx. *Anaesthesia.* 1993;48(12):1075-1078.
27. Liu J, Zhang X, Gong W, et al. Correlations between controlled endotracheal tube cuff pressure and postprocedural complications: a multicenter study. *Anesth Analg.* 2010;111(5):1133-1137.
28. Ryu JH, Han SS, Do SH, Lee JM, Lee SC, Choi ES. Effect of adjusted cuff pressure of endotracheal tube during thyroidectomy on postoperative airway complications: prospective, randomized, and controlled trial. *World J Surg.* 2013;37(4):786-791.
29. Brain AI. The laryngeal mask--a new concept in airway management. *Br J Anaesth.* 1983;55(8):801-805.
30. Wong DT, Tam AD, Mehta V, Raveendran R, Riad W, Chung FF. New supraglottic airway with built-in pressure indicator decreases postoperative pharyngolaryngeal symptoms: a randomized controlled trial. *Can J Anaesth.* 2013;60(12):1197-1203.
31. Townley SA, Herbert PM. Cuff pressure estimation study (abstr. A-1002). *Eur J Anaesthesiol.* 2006;23(Suppl. 37):257-258.
32. Teoh PF, Seet E, Macachor J, Chia N. Accuracy of ProSeal laryngeal mask airway intracuff pressure estimation using finger palpation technique – a prospective, observational study. *Anaesth Intensive Care.* 2012;40(3):467-471.
33. Ong M, Chambers NA, Hullett B, Erb TO, von Ungern-Sternberg BS. Laryngeal mask airway and tracheal tube cuff pressures in children: are clinical endpoints valuable for guiding inflation? *Anaesthesia.* 2008;63(7):738-744.
34. Ghai B, Sethi S, Ram J, Wig J. Cuff filling volumes for pediatric Classic laryngeal mask airways: comparison of clinical end points versus adjusted cuff pressure. *Paediatr Anaesth.* 2013;23(2):122-126.
35. Lumb AB, Wrigley MW. The effect of nitrous oxide on laryngeal mask cuff pressure. In vitro and in vivo studies. *Anaesthesia.* 1992;47(4):320-323.
36. Ouellette RG. The effect of nitrous oxide on laryngeal mask cuff pressure. *AANA J.* 2000;68(5):411-414.
37. Park HS, Han JJ, Kim YJ. The effect of head rotation on efficiency of ventilation and cuff pressure using the PLMA in pediatric patients. *Korean J Anesthesiol.* 2011;61(3):220-224.
38. Miyashiro RM, Yamamoto LG. Endotracheal tube and laryngeal mask airway cuff pressures can exceed critical values during ascent to higher altitude. *Pediatr Emerg Care.* 2011;27(5):367-370.
39. McNamara D. Intracuff monitoring could prevent overinflation injuries. *Anesthesiology News.* 2013;39(4).
40. Slocum AH, Jr., Duffley SC, Gamazo JM, Watral A, Spiegel JE, Slocum AH. Design of a pressure measuring syringe. Abstract presented at the Design of Medical Devices Conference, Minneapolis, MN, USA, 13-15 April 2010. http://web.mit.edu/2.75/projects/DMD_2010_Slocum.pdf. Accessed 15 February, 2016.
41. Maboudi A, Abtahi H, Hosseini M, Tamadon A, Safavi E. Accuracy of endotracheal tube cuff pressure adjustment by fingertip palpation after training of intensive care unit nurses. *Iran Red Crescent Med J.* 2013;15(5):381-384.
42. Lawler PG, Rayner RR. The limitations of the Shiley Pressure Relief Adaptor. *Anaesthesia.* 1982;37(8):865.
43. Spiegel J, Slocum A, Jr., Slocum A, Sr. Measuring endotracheal-tube cuff pressure: an integrated, in-line pressure-measuring syringe (abstr. A1396). Abstract presented at the Annual Meeting of the American Society of Anesthesiologists, San Diego, CA, USA, 16-20 October 2010. <http://www.asaabstracts.com/strands/asaabstracts/abstract.htm?sessionid=C227BB76190CC92BA935491F9981632F?year=2010&index=8&absnum=943>. Accessed 15 February, 2016.
44. Pajjyar AK, Wen Z, Wang H, Ma L, Miao L, Wang G. Comparisons of clinical performance of Guardian laryngeal mask with laryngeal mask airway ProSeal. *BMC Anesthesiol.* 2015;15:69.

[†] The cuff pressure indicator used in the studies cited here was Cuff Pilot Technology, although it was referred to using different brand names.

Teleflex is a global provider of medical technologies designed to improve the health and quality of people's lives.

We apply purpose-driven innovation – a relentless pursuit of identifying unmet clinical needs – to benefit patients and healthcare providers. Our portfolio is diverse, with solutions in the fields of vascular and interventional access, surgical, anesthesia, cardiac care, urology, emergency medicine and respiratory care. Teleflex employees worldwide are united in the understanding that what we do every day makes a difference. For more information, please visit teleflex.com.

Teleflex is the home of Arrow®, Deknatel®, Hudson RCI®, LMA®, Pilling®, Rüschi® and Weck® – trusted brands united by a common sense of purpose.

Corporate Office

Phone +1 610 225 6800, 550 E. Swedesford Road, Suite 400, Wayne, PA 19087, USA

Regional Offices

United States: Phone +1 919 544 8000, Toll Free 866 246 6990, cs@teleflex.com, 3015 Carrington Mill Boulevard, Morrisville, NC 27560, USA

Latin America: Phone +1 919 433 4999, la.cs@teleflex.com, 3015 Carrington Mill Boulevard, Morrisville, NC 27560, USA

International: Phone +353 (0)9 06 46 08 00, orders.intl@teleflex.com, Teleflex Medical Europe Ltd., IDA Business and Technology Park, Dublin Road, Athlone, Co Westmeath, Ireland

Australia/New Zealand 1300 360 226

Austria +43 (0)1 402 47 72

Belgium +32 (0)2 333 24 60

Canada +1 (0) 905 943 9000

China (Shanghai) +86 (0)21 6163 0965

China (Beijing) +86 (0)10 6418 5699

Czech Republic +420 (0)495 759 111

France +33 (0)5 62 18 79 40

Germany +49 (0)7151 406 0

Greece +30 210 67 77 717

India +91 (0)44 2836 5040

Italy +39 0362 58 911

Japan +81 (0)3 6632 3600

Korea +82 2 536 7550

Mexico +52 55 5002 3500

Netherlands +31 (0)88 00 215 00

Portugal +351 22 541 90 85

Singapore (SEA non-direct sales countries) +65 6439 3000

Slovak Republic +421 (0)3377 254 28

South Africa +27 (0)11 807 4887

Spain +34 918 300 451

Switzerland +41 (0)31 818 40 90

United Kingdom +44 (0)1494 53 27 61

For more information, please visit teleflex.com.

Teleflex, the Teleflex logo, Cuff Pilot, LMA, LMA Classic, LMA Fastrach, LMA Flexible, LMA ProSeal, LMA Supreme and LMA Unique are trademarks or registered trademarks of Teleflex Incorporated or its affiliates, in the U.S. and/or other countries. All other trademarks are property of their respective owners.

Information in this document is not a substitute for the product Instructions for Use. The products in this document may not be available in all countries. Please contact your local representative. All data current at time of printing (09/2016). Subject to technical changes without further notice.

© 2016 Teleflex Incorporated. All rights reserved.