

Arrowg+ard Blue Advance™ Midline


Extended dwell*, extended catheter protection

The one midline that delivers extended dwell* and extended catheter protection

When the situation calls for a midline, reach for the pressure-injectable Arrowgard Blue Advance™ Midline. Using an insertion approach you are accustomed to, it offers critical innovations that you and your patients will appreciate.

Arrowgard Blue Advance™ Protection reduces the risk of catheter-related complications— including occlusion, phlebitis, and thrombus accumulation on catheter surfaces— for at least 29 days, giving you the flexibility for extended dwell*.^{2,5,6} It's the midline designed to actually inhibit thrombus formation on catheter surfaces.² Plus, unique features include a GlideThru® Sheath and a Blue FlexTip® feature that helps to minimize vessel damage.¹

Designed to address typical complaints of today's midlines, the Arrowgard Blue Advance™ Midline is the choice that patients, clinicians and hospitals will love for so many reasons.



Protects against thrombus accumulation on catheter surfaces for at least 29 days²



Protects against microbial surface colonization for at least 29 days³




Benefits you—and your patients—will love

Pressure injection
Up to 5 mL/second
through the distal lumen



Single and dual lumen design

Optimal length
15 cm to fit most needs;
cm markings for trimming
to shorter lengths



Unique GlideThru® Sheath

Internal and external protection
Antithrombogenic and antimicrobial
Arrow+ard Blue Advance™
Protection^{2,3}

Staggered exit ports
Reduces risk of mixing
incompatible drugs and
solutions that may create
precipitate⁴




Tapered end for easy advancement
Proprietary Blue FlexTip® Design
minimizes risk of vessel damage¹

Proven protective technology

This innovative midline uses Arrowg+ard Blue Advance™ Protection—proven to provide 29 days of antithrombogenic and antimicrobial catheter protection.^{2,3} Using an initial burst plus sustained release of chlorhexidine, it protects both the internal and external catheter surfaces to reduce these typical catheter-related complications.




**Antimicrobial
Antithrombogenic**




Occlusion

Works continuously to prevent build-up of thrombus, reducing the risk of intraluminal thrombotic occlusion.⁵



Catheter-related Thrombus

Reduces thrombus accumulation on catheter surfaces.²



Phlebitis & Intimal Hyperplasia

An average of 72% less intimal hyperplasia after 29 days and reduction in phlebitis.⁶

Arrowg+ard Blue Advance™ Catheter vs. standard options

In an intravascular *in vivo* model, Arrowg+ard Blue Advance™ Protection demonstrated a total of 92% reduction in thrombus accumulation on the catheter surface when challenged with *Staph aureus* as compared to an uncoated PICC control.² Another *in vivo* model that compared different catheter technologies showed the Arrowg+ard Blue Advance™ catheter to have the least amount of thrombus after indwelling for 30 days.⁶



Unprotected Control Catheter—Significant thrombus and fibrin formation



Bioflo® Catheter with Endexo® Technology—Significant thrombus formation



Arrowg+ard Blue Advance™ Catheter—Minimal thrombus formation²



Your choice matters

As a vascular access specialist, selecting a device that meets your patients' needs—for their specific condition and treatment path—is critical for optimal insertion. That's why you need a range of options for reliable peripheral access.

The more you use midlines, the more you will appreciate the design of the Arrowg+ard Blue Advance™ Midline. It's one more example of how Teleflex is delivering the Right Line For The Right Patient At The Right Time™.



Ordering information

Available in an ergonomically optimized Arrow® ErgoPack® System, with all components needed and arranged in order of use for maximal barriers precautions and clinician safety.

Arrowg+ard Blue Advance™ Midline Basic Catheter Tray • Modified Seldinger

PRODUCT NUMBER	DESCRIPTION	QTY/CASE
PR-41541-BAS	4.5 Fr. x 15 cm Pressure Injectable Arrow® Midline	5
PR-41552-BAS	5.5 Fr. x 15 cm Pressure Injectable Arrow® Midline	5

For more information contact your Teleflex sales representative or call **800.523.8446**.



Rx only

Contraindications: The Arrow+ard Blue Advance™ Midline is contraindicated for patients with known hypersensitivity to chlorhexidine.

*Indwell time is up to 30 days.

References:

1. Rosenbauer KA, Herzer JA. Surface morphology and tensile force at breaking point of different kinds of intravenous catheters before and after usage. *Scan Electron Microsc.* 1981;(Pt 3):125-30.
2. As compared to uncoated catheters, intravascular ovine model inoculated with *Staph aureus*. No correlation between *in vitro/in vivo* testing methods and clinical outcomes have currently been ascertained.
3. *In vitro* data on file 2010. No correlation between *in vitro/in vivo* testing methods and clinical outcomes have currently been ascertained.
4. Collins JL, Lutz RJ. *In vitro* Study of Simultaneous Infusion of Incompatible Drugs in Multilumen Catheters. *Heart & Lung.* 1991; 20(3):271-7.
5. As compared to uncoated catheters, *in vitro* model measuring flush pressure post exposure to human blood. No correlation between *in vitro/in vivo* testing methods and clinical outcomes have currently been ascertained.
6. As compared to uncoated catheters, intravascular ovine model. No correlation between *in vitro/in vivo* testing methods and clinical outcomes have currently been ascertained.

Teleflex is a global provider of medical technologies designed to improve the health and quality of people's lives.

We apply purpose driven innovation – a relentless pursuit of identifying unmet clinical needs – to benefit patients and healthcare providers. Our portfolio is diverse, with solutions in the fields of vascular and interventional access, surgical, anesthesia, cardiac care, urology, emergency medicine and respiratory care. Teleflex employees worldwide are united in the understanding that what we do every day makes a difference. For more information, please visit teleflex.com.

Teleflex is the home of Arrow®, Deknatel®, Hudson RCI®, LMA®, Pilling®, Rüsçh® and Weck® – trusted brands united by a common sense of purpose.

Corporate Office

Phone +1 610 225 6800, 550 E. Swedesford Road, Suite 400, Wayne, PA 19087, USA

Regional Offices

United States: Phone +1 919 544 8000, Toll Free 866 246 6990, cs@teleflex.com, 3015 Carrington Mill Boulevard, Morrisville, NC 27560, USA

Latin America: Phone +1 919 433 4999, la.cs@teleflex.com, 3015 Carrington Mill Boulevard, Morrisville, NC 27560, USA

International: Phone +353 (0)9 06 46 08 00, orders.intl@teleflex.com, Teleflex Medical Europe Ltd., IDA Business and Technology Park, Dublin Road, Athlone, Co Westmeath, Ireland

Australia +1300 360 226

Austria +43 (0)1 402 47 72

Belgium +32 (0)2 333 24 60

Canada +905 943 9004

China (Shanghai) +86 (0)21 6163 0965

China (Beijing) +86 (0)10 6418 5699

Czech Republic +420 (0)495 759 111

France +33 (0)5 62 18 79 40

Germany +49 (0)7151 406 0

Greece +30 210 67 77 717

India +91 (0)44 2836 5040

Italy +39 0362 58 911

Japan +81 (0)3 6632 3600

Korea +82 2 536 7550

Mexico +52 55 5002 3500

Netherlands +31 (0)88 00 215 00

New Zealand +0800 601 100

Portugal +351 22 541 90 85

Singapore (SEA non-direct sales countries) +65 6439 3000

Slovak Republic +421 (0)3377 254 28

South Africa +27 (0)11 807 4887

Spain +34 918 300 451

Switzerland +41 (0)31 818 40 90

United Kingdom +44 (0)1494 53 27 61

For more information, please visit teleflex.com.

Federal Law (USA) restricts these devices to sale by or on the order of a physician.

BioFlo is a trademark and/or registered trademark of AngioDynamics, Inc. Endexo is a trademark of Interface Biologics.

No correlation between *in vitro/in vivo* testing methods and clinical outcomes have currently been ascertained.

Teleflex, the Teleflex logo, Arrow, Arrow+ard Blue Advance, Blue FlexTip, Deknatel, ErgoPack, GlideThru, Hudson RCI, LMA, Pilling, Right Line, Right Patient, Right Time, Rüsçh and Weck are trademarks or registered trademarks of Teleflex Incorporated or its affiliates, in the U.S. and/or other countries. All other trademarks are trademarks of their respective owners.

Information in this material is not a substitute for the product Instructions for Use. Refer to the applicable Instructions for Use for the indications approved in your geography. Please check your local regulatory approval status. Not all products may be available in all countries. Please contact your local representative. Revised: 05/2021.

© 2021 Teleflex Incorporated. All rights reserved. MC-003834 Rev 0.1 LA EN