

Sympro Elite™ Snare

Capture with confidence

Preassembled, one-piece design for rapid deployment


The Sympro Elite Snare eliminates the time needed to assemble separate components, allowing rapid deployment through any 0.035" compatible lumen. When a procedural need arises, the Sympro Elite Snare is ready for immediate use.

Advanced features for added strength and convenience

- Simple 90° loop design remains coaxial to the vessel lumen
- Durable nitinol construction retains its shape and enhances pull strength
- 1:1 torque response for precise positioning and maneuverability
- Radiopaque gold-plated tungsten coil and tip for improved visibility
- Unique locking handle facilitates secure capture


Available in 5, 10, 15, 25, and 35 mm loop diameters for clinical versatility

Redesigned with enhanced pull strength



Sympro Elite™ Snare

MODEL	LOOP DIAMETER	LENGTH	O.D.	BOX QUANTITY
8420	5 mm	150 cm	0.035"	1
8421	10 mm	150 cm	0.035"	1
8422	15 mm	150 cm	0.035"	1
8423	25 mm	150 cm	0.035"	1
8424	35 mm	150 cm	0.035"	1



Please see the Instructions for Use for a complete listing of the indications, contraindications, warnings and precautions.
 CAUTION: Federal Law (U.S.A.) restricts this device to sale by or on the order of a physician.
 Teleflex, the Teleflex logo, and Sympro Elite are trademarks or registered trademarks of Teleflex Incorporated or its affiliates, in the U.S. and/or other countries. Information in this document is not a substitute for the product Instructions for Use. Not all products may be available in all countries. Please contact your local representative.
 Revised: 01/2018. ©All rights reserved. MC-003994 Rev 0