

RÜSCH™

Visualization you can depend on

Rüsch™ Single-Use Fiber-Optic Laryngoscopes

Rüsch™ Polaris™ Blades
Rüsch™ DispoLED™ Handles



Teleflex™
Empowering the future of healthcare

Bright focused light when you need it most.

The Rüscher™ Polaris™ Single-Use Blades with the Rusch™ DispoLED™ Single-Use Handles deliver focused fiber-optic light transmission to aid visualization.




Low-Profile Blade
Designed to help reduce risk of dental trauma

Simple "Click" Locking Mechanism
Provides audible confirmation the blade and compatible handle are secure

"ISO-Green" System
Compatible with ISO 7376 laryngoscope blades and handles

Battery Compartment
All handle components can be recycled after disposal and AAA batteries may be repurposed



Metal Blade Construction
Designed to withstand the forces of intubation with the weight of a reusable blade

Slim Design
for pediatric intubations

Stainless Steel Pin
Built to withstand forces of intubation

Textured Grip
Ergonomically designed finger grips for user comfort and ease of use




A reliable cost efficient single-use solution.

Teleflex offers a broad range of single-use fiber-optic blades and handles to offer a tailored solution to suit various clinician and patient needs.

19 Different Blade Sizes



Macintosh Blades



Miller Blades



Wis-Hipple & Phillips Blades

Handle Options



Medium Handle

Penlight Handle

up to
86%
of reusable handles considered ready-to-use are contaminated with microorganisms despite low-level disinfection¹

up to
82%
of reusable blades considered ready-to-use test positive for residual proteins²



All single-use laryngoscope blade and handle components can be recycled

Ordering Information

Rüsch™ Polaris™ Single-Use Fiber-Optic Laryngoscope Blades

ITEM NUMBER	DESCRIPTION	OVERALL LENGTH (MM)	DISTAL WIDTH (MM)	COLOR CODE	QTY/CASE
Macintosh Blades					
4150100	Mac 0	85.0	10.7	Peach	10
4150110	Mac 1	94.5	9.9	White	10
4150120	Mac 2	111.0	13.0	Blue	10
4150130	Mac 3	133.0	14.0	Yellow	10
4150135	Mac 3.5	145.0	14.5	Red	10
4150140	Mac 4	149.0	13.0	Pink	10
4150150	Mac 5	159.0	14.5	Dark Yellow	10

Miller Blades

4150003	Miller 000	60.0	7.0	Light Green	10
4150002	Miller 00	66.5	12.8	Beige	10
4150001	Miller 0	78.0	12.8	Violet	10
4150010	Miller 1	103.0	12.8	Orange	10
4150015	Miller 1.5	124.0	12.8	Dark Blue	10
4150020	Miller 2	154.0	13.5	Grey	10
4150030	Miller 3	195.0	13.5	Green	10
4150040	Miller 4	205.0	19.0	Lavender	10

Wis-Hipple Blades

4151000	Wis-Hipple 1	102.2	9.2	Light Yellow	10
4151500	Wis-Hipple 1.5	117.8	10.8	Light Grey	10

Phillips Blades

4150011	Phillips 1	113.1	11.6	Brown	10
4150022	Phillips 2	165.1	16.0	Dark Pink	10

Rüsch™ DispoLED™ Single-Use Fiber-Optic Laryngoscope Handles

ITEM NUMBER	DESCRIPTION	CIRCUMFERENCE	DIAMETER	LENGTH	WEIGHT	QTY/CASE
77700	Medium	9.42 cm	2.8 cm	12 cm	73 g	20
22200	Penlight	6.4 cm	2.05 cm	13.5 cm	49 g	20

teleflex.com/laryngoscopes



REFERENCES

- Williams D, Dingley J, Jones C, Berry N. Contamination of laryngoscope handles. *J Hosp Infect.* 2010;74(2):123–128.
- Miller DM, Youkhana I, Karunaratne WU, Pearce A. *Anaesthesia.* 2001;56(11):1069–1072.

INTENDED USE: The laryngoscope is a device intended to allow direct visualization of the laryngeal structure and is used during tracheal intubation.

Rx Only.

Refer to package insert for complete warnings, indications, contraindications, precautions, potential complications and instructions for use.

Teleflex, the Teleflex logo, DispoLED, Polaris, and Rüsch are trademarks or registered trademarks of Teleflex Incorporated and its affiliates, in the U.S. and/or other countries.

© 2025 Teleflex Incorporated. All rights reserved. MC-010783 06 • 25

Teleflex: 3015 Carrington Mill Boulevard, Morrisville, NC 27560
Toll Free: 866 246 6990 • Phone: +1 919 544 8000

Teleflex.com

Teleflex™
Empowering the future of healthcare