

Arrow®

Intra-Aortic 50 cc Balloon Catheters

With New Sizing Guidelines

Meeting your patients' counterpulsation needs

Sizing Guidelines*

Arrow® IAB Catheters	30 cc	40 cc	50 cc
Height (feet/inches)	4' 10" - 5' 4"	5' 4" - 6'	> 5' 8"
Height (centimeters)	147 - 162 cm	162 - 182 cm	> 173 cm
Body Surface Area	< 1.8 m ²	> 1.8 m ²	> 1.8 m ²

*Refer to IFU

50 cc FiberOptix® IAB Catheters

- FiberOptix® Sensor delivers high fidelity AP signal and enables proprietary IntraBeat™ Timing when connected to AutoCAT2WAVE® and AC3 Optimus™ IABP Systems
- 0.027" thin-walled alloy central lumen with 1.5" flex tip design
- Coil-reinforced central lumen to help reduce potential for catheter whip
- Robust polycarbonate bifurcate with large bore connector ensures high rate helium exchange up to 180 bpm

50 cc UltraFlex™ IAB Catheters

- The only coil reinforced 50 cc IAB catheters currently available
- Designed with features:
 - Navigating tortuous anatomy
 - Flexibility
 - Kink resistance to promote optimal gas exchange
- 0.027" thin-walled alloy central lumen with 1.5" flex tip design
- Coil reinforced catheter to help reduce potential for catheter whip
- Robust polycarbonate bifurcate with large bore connector ensures high rate helium exchange up to 180 bpm

50 cc RediGuard® IAB Catheters

- Flexible, soft catheter with stainless-steel braid central lumen
- Centimeter markings on catheter body to facilitate placement

Ordering Information

50 cc FiberOptix® IAB Catheters

PRODUCT NUMBER	CATHETER SIZE	BALLOON VOLUME	INSERTABLE LENGTH	SHEATH LENGTH	CENTRAL LUMIN I.D.	MAXIMUM GUIDEWIRE	CATHETER MATERIAL	BALLOON MEMBRANE LENGTH	INFLATED DIAMETER
IAB-05850-LWS	8 Fr.	50 cc	26.0" (66.0 cm)	6" (15 cm)	0.027"	0.025"	Polyurethane/ Pebax	10.3" (262 mm)	17.0 mm

50 cc UltraFlex™ IAB Catheters

PRODUCT NUMBER	CATHETER SIZE	BALLOON VOLUME	INSERTABLE LENGTH	SHEATH LENGTH	CENTRAL LUMIN I.D.	MAXIMUM GUIDEWIRE	CATHETER MATERIAL	BALLOON MEMBRANE LENGTH	INFLATED DIAMETER
IAB-06850-U	8 Fr.	50 cc	26.0" (66.0 cm)	6" (15 cm)	0.027"	0.025"	Polyurethane/ Stainless Steel	10.3" (262 mm)	17.0 mm

50 cc RediGuard® IAB Catheters

PRODUCT NUMBER	CATHETER SIZE	BALLOON VOLUME	INSERTABLE LENGTH	SHEATH LENGTH	CENTRAL LUMIN I.D.	MAXIMUM GUIDEWIRE	CATHETER MATERIAL	BALLOON MEMBRANE LENGTH	INFLATED DIAMETER
IAB-S850C	8 Fr.	50 cc	26.0" (66.0 cm)	6" (15 cm)	0.027"	0.025"	Polyurethane	10.3" (262 mm)	17.0 mm

Note: All Arrow® IAB catheters are not made with natural rubber latex

Replacement Insertion Kits

PRODUCT NUMBER	DESCRIPTION
IAK-07848	Replacement insertion kit for use with all of Arrow® 50 cc IAB catheters
IAK-02693	Driveline tubing with pre-attached Arrow® 50 cc pump connector for use with Arrow® pump consoles and all Arrow® 50 cc IAB catheters
IAK-02271	Driveline tubing for use with Datascope® Balloon Pumps and Arrow® 50 cc IAB catheters

Additional Accessories

PRODUCT NUMBER	DESCRIPTION
IAA-03270	0.025 in x 175cm PTFE Coated - 3mm "J" Extra Stiff Wire Guide
IB0200	StatLock® Stabilization Device

ACCLEANT

Please see the Instructions for Use for a complete listing of the indications, contraindications, warnings and precautions.

CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.

Teleflex, the Teleflex logo, AC3 Optimus, Arrow, AutoCAT2WAVE, FiberOptix, RediGuard, ProActive Counterpulsation, and UltraFlex are trademarks or registered trademarks of Teleflex Incorporated or its affiliates, in the U.S. and/or other countries. Information in this material is not a substitute for the product Instructions for Use. Not all products may be available in all countries.

All other trademarks are trademarks of their respective owners.

Please contact your local representative. Revised: 12/2020. © 2020 Teleflex Incorporated. All rights reserved. MC-004833 · Rev 1