

DEKLENE® MAXX™ POLYPROPYLENE SUTURE


Strength, give and elasticity

FEEL THE DIFFERENCE WITH DEKLENE MAXX

Knot Percent Elongation at Break

The stress that causes suture breakage occurs during knotting, which makes the comparative evaluation for Knot Percent Elongation at Break most relevant when assessing suture performance.

DEKLENE MAXX VS. LEADING COMPETITOR*



*Test reports on file with Teleflex.

DEKLENE MAXX


Maximized elongation for warning before breakage



INCREASED KNOT SECURITY:

THE OUTPERFORMING POLYPROPYLENE SUTURE

- Enhanced elongation provides tactile feedback before breakage
- Unmatched suture strength
- Superior knot security
- Reduced memory



PREMIUM XT NEEDLES:

THE ANASTOMOTIC PROBLEM-SOLVER

XT needle series enable easy penetration of calcified vessels

- 300 series stainless steel for superior strength and ductility
- Square body design enhances strength and maintains stability in the needle holder



PRE-PLEDGETED:

REDUCE SET-UP TIME

- Increase efficiency
- Reduce occurrence of frayed pledgets



SINGLE STRAND, STRAIGHT PACK:

PACKAGING TO MEET YOUR NEEDS

- Single strand straight packs further reduce suture memory
- Convenient direct-arm access

TO PLACE AN ORDER CONTACT CUSTOMER SERVICE AT 1.866.246.6990

ORDERING INFORMATION

NEEDLE		LENGTH	SUTURE SIZE					STRANDS/ PACK	PACKS/ BOX	
			6-0	5-0	4-0	3-0	2-0			
T-2 17 mm										
	3/8 Circle	24"			D-737K	D-738K		1	36	
		36"		D-7069K	D-7068K	D-6784K		1	36	
				CD-7068K						
				GD-7068K						
				DM-4570 ^L						
				HD-7068M4K			4	12		
	48"		D-6814K	D-6781K	D-6783K		1	36		
T-10 12.7 mm										
	1/2 Circle	18"	D-7375K	D-7374K				1	36	
		24"	D-7275K	D-7274K				1	36	
		30"	D-7175K	D-7174K				1	36	
			D-7175M4K	D-7174M4K				4	12	
		36"		D-8205K				1	36	
T-3 16.2 mm										
	1/2 Circle	36"			D-6478K ^C			1	36	
AT-2 17.4 mm										
	1/2 Circle	18"		D-7064K	D-6774K			1	36	
		24"		D-6150K				1	36	
		36"		D-6151K	D-7065K	D-7066K		1	36	
					CD-7065K	CD-7066K				
					GD-7065K					
					DM-6749 ^M					
				D-6151M4K	D-7065M4K	D-7066M4K		4	12	
					BD-7065M4K	CD-7066M4K				
				CD-7065M4K						
		48"			D-5793K	D-6813K		1	36	
T-5 19.8 mm										
	1/2 Circle	24"			D-6809K			1	36	
		30"			CD-6149M4K			4	12	
C-13 21.1 mm										
	1/2 Circle	30"				D-7227K		1	36	
		36"			D-7126K	D-7127K		1	36	
					CD-7126K	CD-7127K				
					CD-7126M4K					4
		48"			D-6812K		D-6785K	1	36	
C-2 25.9 mm										
	1/2 Circle	36"			D-7041K	D-7042K	D-7044K	1	36	
					CD-7041K	CD-7042K				
					GD-7041K	GD-7042K				
					DM-6775 ^L	DM-4760 ^J				
					DM-6748 ^M	DM-4567 ^N				
					DM-5038 ^I					
					D-7041M4K	D-7042M4K				4
			CD-7041M4K	CD-7042M4K						
		48"			D-5792K	D-5795K		1	36	
HR36+ 36 mm										
	1/2 Circle	30"				D-6806K		1	36	

NEEDLE		LENGTH	SUTURE SIZE				STRANDS/ PACK	PACKS/ BOX
			8-0	7-0	6-0	5-0		
XT-7	6.5 mm							
	3/8 Circle	18"	D-7501K				1	36
			D-7501M1K				1*	12
			D-7501M2K	D-6063M2K			2	12
			D-7501M4K				4	12
	24"	D-7504K	D-6064K			1	36	
		D-7504M1K				1*	12	
		D-7504M2K	D-6064M2K			2	12	
		D-7504M4K	D-6064M4K			4	12	
	30"	D-6551M1K				1*	12	
			D-6247M4K			4	12	
XT-6.5	8 mm							
	3/8 Circle	18"	D-4656M1K				1*	12
				D-5521K			1	36
			D-4656M2K	D-5521M2K			2	12
				D-5521M4K			4	12
		24"	D-7505K	D-7506K			1	36
			D-7505M1K	D-7506M1K			1*	12
			D-7505M2K	D-7506M2K			2	12
			D-7505M4K	D-7506M4K			4	12
	30"	D-4666K	D-7507K			1	36	
		D-4666M2K	D-7507M2K			2	12	
			D-7507M4K			4	12	
XT-6	9.5 mm							
	3/8 Circle	18"		D-7062K			1	36
				D-7062M4K	D-793M4K		4	12
		24"		D-7004K	D-5007K		1	36
				D-7004M1K			1*	12
				D-7004M2K			2	12
				D-7004M4K	D-5007M4K		4	12
		30"		D-7063K	D-7070K		1	36
				D-7063M1K			1*	12
				D-7063M2K			2	12
				D-7063M4K	D-7070M4K		4	12
		36"		D-6365M4K			4	12
XT-1	13 mm							
	3/8 Circle	18"			D-655K		1	36
		24"			D-7016K	D-7060K	1	36
					D-7016M2K	D-7060M2K	2	12
					D-7016M4K	D-7060M4K	4	12
			30"			D-7076K		1
					D-730K**		1	36
					D-7076M1K		1*	12
					D-7076M2K		2	12
					D-7076M4K		4	12
		36"				D-7017K	D-7061K	1
						D-7061M2K	2	12
						D-7061M4K	4	12
AT-6	10 mm							
	1/2 Circle	24"	D-4523K				1	36
				FD-4522M2K			2	12
		30"	D-4446K				2	36

PREFIXES/FOOTNOTES:

B = SOFT 1/8" X 1/4" C = SOFT 3/8" X 3/16" F = FIRM 1/8" X 1/8" G = FIRM 1/8" X 1/4" H = FIRM 3/8" X 3/16" J = SOFT 1/2" X 1/4" L = SOFT 1/2" X 5/16"
 M = SOFT 1/4" X 3/16" N = SOFT 5/8" X 3/8" *SINGLE STRAIGHT PACK **ONE NEEDLE

Teleflex is a global provider of medical products designed to enable healthcare providers to protect against infections and improve patient and provider safety. The company specializes in products and services for vascular access, respiratory, general and regional anesthesia, cardiac care, urology and surgery. Teleflex also provides specialty products for device manufacturers.

The Teleflex family of brands includes ARROW®, BEERE MEDICAL®, DEKNATEL®, GIBECK®, HUDSON RCI®, KMEDIC®, PILLING®, PLEUR-EVAC®, RÜSCH®, SHERIDAN®, SMD®, TAUT®, TFX OEM®, VASONOVA™ and WECK®, all of which are trademarks or registered trademarks of Teleflex Incorporated or its affiliates. Teleflex Incorporated (NYSE:TFX) has annual revenues of approximately \$1.5 billion in revenue and customers in more than 130 countries.

For detailed information, see teleflex.com

TELEFLEX MEDICAL

PO Box 12600

Research Triangle Park, NC 27709

Toll Free 866.246.6990

Phone 919.544.8000

Intl 919.433.8088

Teleflex, Deklene, Deknatel and Maxx are trademarks or registered trademarks of Teleflex Incorporated or its affiliates. © 2012 Teleflex Incorporated. All rights reserved. 2012-0841