

Pleur-evac™
Chest Drainage Systems



Thoracic Catheter & Accessory Data Sheet



To view Pleur-evac™
Systems and Accessories
go to myTeleflex.com



Thoracic Catheters and Accessories



PE100



PE101



PE102



PE103



PE104-08



PE105

Connectors

REF.	DESCRIPTION	CONNECTOR OUTER DIAMETER (O.D.)	PATIENT TUBING PAIRING	STRAIGHT/RIGHT ANGLE COMPATIBILITY	TROCAR CATHETER COMPATIBILITY	QTY
PE100	3/8" Universal Connector	9.5mm / 14.0mm	Pre-connected to all adult units & S-1120	24-40 Fr.	28-32 Fr.	20
PE101	3/8" Y Connector	9.5 mm	For y-connecting 3/8" (9.5mm) patient tubing; pre-connected to S-1120			20
PE102	5-in-1 Connector	7.0mm / 14.0mm	Individually packaged with infant units A-6020, A-8020, & S-1130; Fits 1/4" & 3/8" tubing	16-40 Fr.	20-32 Fr.	20
PE103	5mm x 10mm Straight Connector	3.5mm / 9.0mm	Pre-connected to infant units A-6020 & A-8020; Fits 1/4" tubing	12-24 Fr.	10-24 Fr.	20
PE104-08	Infant Connector	5.0mm / 6.5mm	Pre-connected to infant unit S-1130; Fits 3/8" tubing	12-20 Fr.	10-16 Fr.	20
PE105	Patient Tubing Plug		For disposal of a Pleur-evac™ Chest Drainage Device or occluding a single drainage site of a dual drainage unit.			20

Additional Tubing

S-1120	Dual Patient Drainage Tube	Compatible with Pleur-evac™ Chest Drainage systems containing quick-disconnect tubing; equipped with two 6ft-lengths of patient tubing (not made with natural latex); joined together by y-connector PE101; pre-connected with two PE100 universal, b connectors	6
--------	----------------------------	------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------	---



Straight Catheters

REF.	DESCRIPTION	QTY
DSTC-12S	Straight – Soft Thoracic Catheter, 18-3/8", 12 Fr.	10
DSTC-16S	Straight – Soft Thoracic Catheter, 18-3/8", 16 Fr.	10
DSTC-20S	Straight – Soft Thoracic Catheter, 18-3/8", 20 Fr.	10
DSTC-24S	Straight – Soft Thoracic Catheter, 18-3/8", 24 Fr.	10
DSTC-28S	Straight – Soft Thoracic Catheter, 18-3/8", 28 Fr.	10
DSTC-32S	Straight – Soft Thoracic Catheter, 18-3/8", 32 Fr.	10
DSTC-36S	Straight – Soft Thoracic Catheter, 18-3/8", 36 Fr.	10
DSTC-40S	Straight – Soft Thoracic Catheter, 18-3/8", 40 Fr.	10

Right Angle Catheters

REF.	DESCRIPTION	QTY
DRAC-16S	Right Angle – Soft Thoracic Catheter, 18-3/8", 16 Fr.	10
DRAC-20S	Right Angle – Soft Thoracic Catheter, 18-3/8", 20 Fr.	10
DRAC-24S	Right Angle – Soft Thoracic Catheter, 18-3/8", 24 Fr.	10
DRAC-28S	Right Angle – Soft Thoracic Catheter, 18-3/8", 28 Fr.	10
DRAC-32S	Right Angle – Soft Thoracic Catheter, 18-3/8", 32 Fr.	10
DRAC-36S	Right Angle – Soft Thoracic Catheter, 18-3/8", 36 Fr.	10
DRAC-40S	Right Angle – Soft Thoracic Catheter, 18-3/8", 40 Fr.	10

Trocar Catheters

REF.	DESCRIPTION	QTY
DTRC-10S	Thoracic Catheter Trocar, 7-1/8", 10 Fr.	10
DTRC-12S	Thoracic Catheter Trocar, 8-1/8", 12 Fr.	10
DTRC-16S	Thoracic Catheter Trocar, 8-1/2", 16 Fr.	10
DTRC-20S	Thoracic Catheter Trocar, 14", 20 Fr.	10
DTRC-24S	Thoracic Catheter Trocar, 14", 24 Fr.	10
DTRC-28S	Thoracic Catheter Trocar, 14", 28 Fr.	10
DTRC-32S	Thoracic Catheter Trocar, 14", 32 Fr.	10