


182911
Extended-Length
Judd Allis Forceps



Pilling™ Extended-Length Judd Allis Forceps

A Versatile Instrument to Support Your Vascular Needs

Now available in the Pilling™ Line of CVT Instrumentation

Features

- Extended length adaptation of existing 7-1/2" Judd Allis Forceps (182905)
 - Previously only available through Custom request process
- Additional 3-3/4" shaft length for deeper reach
- Low profile jaw for atraumatic approach
- Versatility for a wide range of procedures including: Thoracic, Vascular, and even General, GYN, and trauma.

Judd Allis Family

ITEM CODE	DESCRIPTION	OVERALL LENGTH	JAW
182911	NEW Extended-Length Judd Allis Forceps	11" (28 cm)	1-3/4" (4.5 cm), 3x4 teeth with horizontal serrations
182905	Judd Allis Forceps	7-1/4" (19 cm)	1-3/4" (4.5 cm), 3x4 teeth with horizontal serrations
182900	Judd Allis Forceps	5-7/8" (15 cm)	2" (4.9 cm), 3x4 teeth only

Other Specialized Pilling™ Vascular Instruments

Heitmiller Mini Thoracotomy Retractor

This retractor is designed to provide exposure for video-assisted thoracoscopy, open lung biopsies, transaxillary incisions and MIDCAB procedures in which mini-thoracotomy incisions are used. With its small scale it is also utilized to mobilize lung adhesions during video thoracotomies.



341310
Heitmiller Mini
Thoracotomy Retractor

Kaiser-Pilling™ Thoracoscopy Instrumentation

These instruments deliver the versatility you require and the feel, balance, and weight to which you are accustomed. Featuring traditional instrument patterns that have been elongated and shaped to meet the special needs of the surgeon, providing an extended reach.



841032
Kaiser-Pilling™
Ring Forceps
Short Jaw C-Shape

ITEM CODE	DESCRIPTION
841010	Kaiser-Pilling™ Ring Forceps Long Jaw Straight
841011	Kaiser-Pilling™ Ring Forceps Long Curved Jaw
841012	Kaiser-Pilling™ Ring Forceps Long Jaw C-Shape
841013	Kaiser-Pilling™ Ring Forceps Long Jaw Curved Side
841014	Kaiser-Pilling™ Ring Forceps Long Jaw Curved Handle
841030	Kaiser-Pilling™ Ring Forceps Short Jaw Straight
841031	Kaiser-Pilling™ Ring Forceps Short Jaw Curved
841032	Kaiser-Pilling™ Ring Forceps Short Jaw C-Shape
841033	Kaiser-Pilling™ Ring Forceps Short Jaw Curved Side
841050	Kaiser-Pilling™ Ring Forceps Straight Atraugrip

ITEM CODE	DESCRIPTION
841051	Kaiser-Pilling™ Ring Forceps S-Shape Atraugrip
841052	Kaiser-Pilling™ Ring Forceps Curved Jaw Atraugrip
841100	Kaiser-Pilling™ Biopsy Forceps, 4 mm Cup Malleable
841101	Kaiser-Pilling™ Biopsy Forceps, 4 mm Cup Right Curved
841300	Kaiser-Pilling™ Spatula Slight Curve
841301	Kaiser-Pilling™ Spatula Acute Curve
841350	Kaiser-Pilling™ Yankauer With Thumb Valve
841351	Kaiser-Pilling™ Andrews Large Diameter Tube
841400	Kaiser-Pilling™ Hook Blunt



To order the NEW Extended-Length Judd Allis, visit [myTeleflex.com](https://www.myTeleflex.com)

Teleflex, the Teleflex logo, and Pilling are trademarks or registered trademarks of Teleflex Incorporated or its affiliates, in the U.S. and/or other countries. Information in this document is not a substitute for the product Instructions for Use. The products in this document may not be available in all countries. Please contact your local representative. All data current at time of printing (01/2024). Subject to technical changes without further notice.

© 2024 Teleflex Incorporated. All rights reserved. MC-009494