




532221
NEW Pilling™ Dedo
Laser Laryngoscope



532196
NEW Garrett-Ossoff-Pilling™
Laryngoscope



50SC06L01 Pilling™ MicroLaryngoscopy Solid Container System

Pilling™
Microlaryngoscopy
Organization Configurations

Organization Options


Protect and organize your full microlaryngoscopy procedural set in one container

Solid Container



- Available in Microlaryngoscopy configuration only (50SC06L01): 27" x 11" x 6"
- Features three basket levels: bottom level holds six laryngoscopes, top two levels store instrumentation
- Solid bottom is compatible with Steam, EtO, and Sterizone VP4 sterilization methods
- Utilizes round sterilization filters

Wrapped Container



- Microlaryngoscopy configuration (P88097): 23" x 11" x 5"
- Features two tiers: bottom tier holds six laryngoscopes, top tier stores instrumentation
- Perforated container requires wrapping for sterilization
- Also available in Bronchoscopy, Esophagoscopy, and Panendoscopy configurations

Microlaryngoscopy Configurations

Solid Containers

Solid-bottom container system.

Three basket levels hold 6 laryngoscopes, plus instruments and light cables.

<input type="checkbox"/>	50SC06L01	Pilling™ Micro Laryngoscopy Container, 27" x 11" x 6", with Starter Kit Starter Kit includes: 40 disposable filters, 20 tamper evident seals, 20 load cards, and 20 load indicator strips
--------------------------	-----------	--

Wrapped Containers

<input type="checkbox"/>	P88097	Double-Decker ENT Tray w/Inserts and Standard Mats, 23" x 11" x 5"
<input type="checkbox"/>	P88044	Double-Decker ENT Tray w/Inserts Only, 23" x 11" x 5"
<input type="checkbox"/>	087396 x2	Ultrasoft Silicone Spiked Mat, 22" L x 11" W – to be paired with P88044

Laryngoscopes

All scopes include Fiber Optic Light Carriers.

<input type="checkbox"/>	532221	NEW Adult Dedo Laser Micro Laryngoscope, 19 mm x 178 mm
<input type="checkbox"/>	532206	NEW Jako Type 3 Laryngoscope, 29 mm x 178 mm
<input type="checkbox"/>	532196	NEW Garrett-Ossoff-Pilling Laryngoscope, 25 mm x 188 mm
<input type="checkbox"/>	532191	NEW Ossoff-Pilling Laser Laryngoscope, 25 mm x 168 mm
<input type="checkbox"/>	532190	NEW Ossoff-Pilling Laser Laryngoscope, 19 mm x 168 mm
<input type="checkbox"/>	532010	NEW Holinger Anterior Commissure Laryngoscope, 10 mm x 170 mm
<input type="checkbox"/>	532215	NEW Lindholm-style Laser Laryngoscope, 38 mm x 150 mm
<input type="checkbox"/>	532217	NEW Lindholm-style Laser Laryngoscope, 38 mm x 170 mm

Laryngeal Instruments

Cup Forceps – Standard (22 cm – 28 cm)

<input type="checkbox"/>	506455	Jako 2 mm Biopsy Straight, 22 cm
<input type="checkbox"/>	506460	Jako 2 mm Biopsy Angled Left, 22 cm
<input type="checkbox"/>	506461	Jako 2 mm Biopsy Angled Right, 22 cm
<input type="checkbox"/>	506462	Jako 2 mm Angled Up, 22 cm
<input type="checkbox"/>	505100	Jackson 4 mm Straight, 28 cm
<input type="checkbox"/>	505104	Jackson 4 mm Angled Up, 28 cm
<input type="checkbox"/>	505102	Jackson 6 mm Straight, 28 cm
<input type="checkbox"/>	505106	Jackson 6 mm Angled Up, 28 cm

Micro-Cup Forceps (23 cm)

<input type="checkbox"/>	506862	Phono 1 mm Straight
<input type="checkbox"/>	506864	Phono 1 mm Angled Up
<input type="checkbox"/>	P28526	Phono 1 mm Right
<input type="checkbox"/>	P28527	Phono 1 mm Left

Grasping Forceps – Standard (22 cm)

<input type="checkbox"/>	506470	Jako Alligator Cvd Left
<input type="checkbox"/>	506471	Jako Alligator Cvd Right
<input type="checkbox"/>	506472	Jako Alligator Straight
<input type="checkbox"/>	506950	Laryngeal Triangular Jaw Right
<input type="checkbox"/>	506952	Laryngeal Triangular Jaw Left

Micro-Grasping Forceps (23 cm)

<input type="checkbox"/>	506856	Phono 1 mm Double Action with Suction
<input type="checkbox"/>	506858	Phono 1 mm Angled Up
<input type="checkbox"/>	506860	Phono 1 mm Straight
<input type="checkbox"/>	P28531	Phono 1 mm Right
<input type="checkbox"/>	P28532	Phono 1 mm Left

Scissors – Standard (22 cm)

<input type="checkbox"/>	506475	Jako Curved Left
<input type="checkbox"/>	506476	Jako Curved Right
<input type="checkbox"/>	506477	Jako Slight Up Angle
<input type="checkbox"/>	506478	Jako Straight Vertical
<input type="checkbox"/>	506479	Jako Straight Horizontal

Suction Tubes (20 cm – 35 cm)

<input type="checkbox"/>	507010	12 Fr Jackson Open End, 20 cm
<input type="checkbox"/>	507030	12 Fr Jackson Velvet Eye, 20 cm
<input type="checkbox"/>	507065	12 Fr Holinger Velvet Eye, 35 cm
<input type="checkbox"/>	507060	9 Fr Holinger Velvet Eye, 20 cm
<input type="checkbox"/>	507122	Jako 7 Fr with Thumb Control
<input type="checkbox"/>	507121	Jako Suction Tube 6 Fr with Thumb Control
<input type="checkbox"/>	507120	Jako 5 Fr with Thumb Control
<input type="checkbox"/>	507045	Insulated 12 Fr Suction Cautery Tube, 22 cm w/10 ft Bovie Cord & Plug

Knives, Probes, and Elevators (27 cm)

<input type="checkbox"/>	506673	Jako Blunt Probe Angled
<input type="checkbox"/>	506669	Jako Blunt Probe Straight
<input type="checkbox"/>	506880	Phono Sickle Knife Straight
<input type="checkbox"/>	506882	Phono Sickle Knife Curved
<input type="checkbox"/>	506830	Phono Probe Elevator 30° Blunt
<input type="checkbox"/>	506832	Phono Probe Elevator 60° Blunt
<input type="checkbox"/>	506888	Universal Ring Handle with Camlock

Ventilation

<input type="checkbox"/>	502600	Venturi Ventilating Needle, Blunt 14 Gauge, Right Angle, 4"
<input type="checkbox"/>	502594	Venturi Ventilating Needle, Blunt 14 Gauge, Straight, 5-3/8"
<input type="checkbox"/>	522240	Clip for Venturi Ventilating Needles

Accessories

Laryngoscope Holders / Table Mounted Tray

<input type="checkbox"/>	502245	Lewy-style Laryngoscope Holder
<input type="checkbox"/>	502248	Tobey Laryngoscope Holder with Quick Release
<input type="checkbox"/>	502260	Boston University Suspension System
<input type="checkbox"/>	502261	Boston University Laryngoscope Handle Extension
<input type="checkbox"/>	502265	Dedo Table Mounted Mayo Tray


Light Cable

<input type="checkbox"/>	P358PRA	Pilling™ to Pilling™ Blue Fiber Optic Cable, 8 ft
<input type="checkbox"/>	A358PRA	ACMI to Pilling™ Blue Fiber Optic Cable, 8 ft

Other Wrapped Container Configurations

Bronchoscopy

<input type="checkbox"/>	P37331	Single-Tier ENT Tray w/o Brackets, 23" x 11" x 3.5"
<input type="checkbox"/>	P37306	x3 Universal Keyhole Bracket 15 Position, 10"
<input type="checkbox"/>	P37305	x2 Accessory Box w/Hardware, 3.5" x 1.5" x 1"
<input type="checkbox"/>	087396	Ultrasoft Silicone Spiked Mat, 22" L x 11" W



Bronchoscopy configuration in tray P37331, 23" x 11" x 3.5"

Esophagoscopy

<input type="checkbox"/>	P37331	Single-Tier ENT Tray w/o Brackets, 23" x 11" x 3.5"
<input type="checkbox"/>	P37306	x2 Universal Keyhole Bracket 15 Position, 10"
<input type="checkbox"/>	087396	Ultrasoft Silicone Spiked Mat, 22" L x 11" W




Esophagoscopy configuration in tray P37331, 23" x 11" x 3.5"

Esophagoscopy and Bronchoscopy Instrumentation

<input type="checkbox"/>	P37304	Single-Tier ENT Tray w/o Brackets, 27" x 11" x 3.5"
<input type="checkbox"/>	P37306	x2 Universal Keyhole Bracket 15 Position, 10"
<input type="checkbox"/>	087396	Ultrasoft Silicone Spiked Mat, 22" L x 11" W

Extra-Long - For storage of 60cm Esophageal Instruments

<input type="checkbox"/>	P37319	Extra-Long Tray w/Brackets, 31" x 6.5" x 3"
--------------------------	--------	---




Bronchoscopy & Esophagoscopy Instrument configuration in tray P37304, 27" x 11" x 3.5"

Panendoscopy

<input type="checkbox"/>	P37304	Single-Tier ENT Tray w/o Brackets, 27" x 11" x 3.5"
<input type="checkbox"/>	P37338	Tray Shelf, 21.5" x 10.5" x 2.25"
<input type="checkbox"/>	P37545	x4 "L" Brackets to Support Shelf
<input type="checkbox"/>	P37306	x5 Universal Keyhole Bracket 15 Position, 10"
<input type="checkbox"/>	P37379	x2 Keyhole Bracket, Single, 1.5"
<input type="checkbox"/>	P37305	Accessory Box w/Hardware, 3.5" x 1.5" x 1"
<input type="checkbox"/>	087396	x2 Ultrasoft Silicone Spiked Mat, 22" L x 11" W



Extra-Long Esophagoscopy Instrument configuration in tray P37319, 31" x 6.5" x 3"



Panendoscopy configuration in tray P37304, 27" x 11" x 3.5" with tray shelf P37338 (shelf provides a double decker option)

Note:

1. Refer to "Tray & Components" checklist for all components necessary to build pictured configurations
2. Assembly required. Sales Rep assistance available upon request.

Teleflex is a global provider of medical technologies designed to improve the health and quality of people's lives. We apply purpose driven innovation – a relentless pursuit of identifying unmet clinical needs – to benefit patients and healthcare providers. Our portfolio is diverse, with solutions in the fields of vascular access, interventional cardiology and radiology, anesthesia, emergency medicine, surgical, and urology. Teleflex employees worldwide are united in the understanding that what we do every day makes a difference. For more information, please visit teleflex.com.

Teleflex is the home of Arrow™, Deknatel™, LMA™, Pilling™, QuikClot™, Rüschi™, UroLift™, and Weck™ – trusted brands united by a common sense of purpose.

Corporate Office

Phone +1 610 225 6800, 550 E. Swedesford Road, Suite 400, Wayne, PA 19087, USA

Regional Offices

United States: Phone +1 919 544 8000, Toll Free 866 246 6990, cs@teleflex.com, 3015 Carrington Mill Boulevard, Morrisville, NC 27560, USA

Latin America: Phone +1 919 433 4999, la.cs@teleflex.com, 3015 Carrington Mill Boulevard, Morrisville, NC 27560, USA

International: Phone +353 (0)9 06 46 08 00, orders.intl@teleflex.com, Teleflex Medical Europe Ltd., IDA Business and Technology Park, Dublin Road, Athlone, Co Westmeath, Ireland

Australia 1300 360 226

Austria +43 (0)1 402 47 72

Belgium +32 (0)2 333 24 60

Canada +1 (0)800 387 9699

China (Shanghai) +86 (0)21 6163 0965

China (Beijing) +86 (0)10 6418 5699

Czech Republic +420 (0)495 759 111

France +33 (0)5 62 18 79 40

Germany +49 (0)7151 406 0

Greece +30 210 67 77 717

India +91 (0)44 2836 5040

Italy ++39 0362 58911

Japan +81 (0)3 6632 3600

Korea +82 2 536 7550

Mexico +52 55 5002 3500

Netherlands +31 (0)88 00 215 00

New Zealand 0800 601 100

Portugal +351 22 541 90 85

Singapore (SEA non-direct sales countries) +65 6439 3000

Slovak Republic +421 (0)3377 254 28

South Africa +27 (0)11 807 4887

Spain +34 918 300 451

Switzerland +41 (0)31 818 40 90

United Kingdom +44 (0)1494 53 27 61

For more information, please visit teleflex.com.

Teleflex, the Teleflex logo, Arrow, Deknatel, LMA, Pilling, QuikClot, Rüschi, UroLift, and Weck are trademarks or registered trademarks of Teleflex Incorporated or its affiliates in the U.S. and/or other countries.

Information in this material is not a substitute for the product Instructions for Use. Not all products may be available in all countries. Please contact your local representative. Revised: 04/2024.

© 2024 Teleflex Incorporated. All rights reserved.

MC-007620 Rev 1