



Welcome to the Family!

Pilling®  
**Lindholm-Style Laser Laryngoscope**  
 The Newest Addition to the Pilling® Family


**Lindholm-Style Laser Laryngoscope**

The new Lindholm-style laser laryngoscope ushers in a new era of Pilling® laryngoscopes. Each Lindholm is compatible with a single light carrier, smoke evacuation tube, and jet ventilation needle. All accessories come in a modular design allowing the surgeon to tailor the laryngoscope to the patient and condition being treated.

**Features and Benefits**

- Modular accessories
- Color coded accessories
- Laser compatible
- Jet ventilation compatible
- Visualization of posterior commissure, epiglottis, and cricoid region

# Pilling® Lindholm-Style Laser Laryngoscope




## Pilling® Lindholm-Style Laryngoscopes

PRODUCT CODE	DESCRIPTION	SIZE	MINIMAL DISTAL OPENING	MOUTH PROXIMAL DIMENSIONS
532215	Adult 15 cm Full System	150 mm	18 mm	38 mm x 24 mm
532215S	Scope Only	150 mm		
532217	Adult 17 cm Full System	170 mm		
532217S	Scope Only	170 mm		

## Individual or Replacement Accessories

PRODUCT CODE	DESCRIPTION
537215	Light Carrier
538215	Smoke Evac Tube
539215	Jet Ventilation Needle
Color Code	Blue



To view Pilling® Laryngoscope products go to:  
[myTeleflex.com](https://myTeleflex.com)

Teleflex, the Teleflex logo, and Pilling are trademarks or registered trademarks of Teleflex Incorporated or its affiliates in the U.S. and/or other countries. Information in this material is not a substitute for a product's Instructions for Use. Not all products may be available in all countries. Please contact your local representative. Revised: 04/2023  
© 2023 Teleflex Incorporated. All rights reserved. MC-008869 Rev 0