


The quality of care you demand for your patients

Arrow® FlexTip Plus® Epidural Catheter now with NRFit™ connectors.

The Arrow® FlexTip Plus® Epidural Catheter from Teleflex is clinically proven to significantly reduce complications in the epidural space.¹⁻⁶ This unique polyurethane catheter is designed to provide the precise balance of characteristics needed to get ease of placement, reliable patency, and optimum indwelling characteristics. This technology is available in both open-tip and closed-tip catheter configurations, so you can achieve the drug dispersion you prefer with the quality of care you demand for your patients.



Teleflex has adopted ISO 80369-6 standards to support efforts to improve patient safety in neuraxial settings. The NRFit™ connector is designed to reduce the possibility of misconnection between devices that are not intended to connect.

Reliable patency

The coil reinforcement makes the FlexTip Plus® Catheter kink resistant to better maintain patency, even in some situations that might cause ordinary nylon catheters to fail.^{2,3}

Less inadvertent penetration

The soft, flexible tip is designed to deflect when it contacts the dura or blood vessels, which significantly reduces the risk of vessel puncture¹⁻⁵ and paresthesia.^{2,6}

The results are enhanced comfort and safety for the patient and greater ease of use and reliability for the anesthesia provider.

Arrow® Epidural Catheterization Kits and Accessories with NRFit™ connectors

19 Ga. Arrow® FlexTip Plus® Open-Tip Epidural Catheter Trays

ITEM NUMBER	CONFIGURATION	EPIDURAL NEEDLE	LOSS OF RESISTANCE SYRINGE	MEDICATIONS	ADVANCER TUBE	UNITS/ BOX	COMPARABLE PRODUCT WITH STANDARD CONNECTOR
OT19TAGSK	Kit	17 Ga. x 3-1/2" Tuohy, Winged	5ml, Slip, Glass	•		10	AK-05502
OT19TAGSKA	Kit	17 Ga. x 3-1/2" Tuohy, Winged	5ml, Slip, Glass	•	•	10	SJ-05501
OT19TAPSK	Kit	17 Ga. x 3-1/2" Tuohy, Winged	10ml, Slip, Plastic	•		10	AK-05501
OT19TCGLK	Kit	17 Ga. x 3-3/8" Tuohy, Removable Wings	5ml, Lock, Glass	•		10	AK-05503-L
OT19HAGSKA	Kit	17 Ga. x 3-7/8" Hustead, Removable Wings	5ml, Slip, Glass	•	•	10	WJ-05401
OT19TKPSS	Set	17 Ga. x 3-1/2" Tuohy, Winged	10ml, Slip, Plastic			20	JC-05400-B
OT19HBPSS	Set	17 Ga. x 3-7/8" Hustead, Removable Wings	10ml, Slip, Plastic			20	JC-05400-E

19 Ga. Arrow® FlexTip Plus® Open-Tip Combined Spinal/Epidural Catheter Trays

ITEM NUMBER	CONFIGURATION	EPIDURAL NEEDLE	SPINAL NEEDLE	LOSS OF RESISTANCE SYRINGE	MEDICATIONS	ADVANCER TUBE	UNITS/ BOX	COMPARABLE PRODUCT WITH STANDARD CONNECTOR
CSE19TAGSKG	Kit	17 Ga. x 3-1/2" Tuohy, Winged	26 Ga. x 5" Gertie Marx® Spinal	5ml, Slip, Glass	•		10	AK-05560

19 Ga. Arrow® TheraCath® Epidural Catheter Trays

ITEM NUMBER	CONFIGURATION	EPIDURAL NEEDLE	LOSS OF RESISTANCE SYRINGE	MEDICATIONS	ADVANCER TUBE	UNITS/ BOX	COMPARABLE PRODUCT WITH STANDARD CONNECTOR
TC19HAGLN	Intermediate Kit	17 Ga. x 3-7/8" Hustead	5ml, Slip, Glass			10	AK-05000

Catheters and Accessories with NRFit™ Connectors

ITEM NUMBER	CONFIGURATION	UNITS/ BOX	COMPARABLE PRODUCT WITH STANDARD CONNECTOR
ECTC19	19 Ga. TheraCath® Epidural Catheter with Stylet, SnapLock™ Adapter, and Threading Assist Device	10	EC-05000
SL19N	SnapLock™ Adapter with NRFit™ connector	50	AS-05500

For more information
Visit flextipplus.com

19 Ga. FlexTip Plus® Open-Tip Epidural Catheterization Kit with NRFit™ connectors Includes:

- FlexTip Plus® Epidural Catheter with Thread Assist Device with NRFit™ Connector
- Epidural Needle with NRFit™ Connector
- Loss of Resistance Syringe with NRFit™ Connector
- SnapLock™ Catheter Syringe Adapter with NRFit™ Connector
- Injection Needle: 25 Ga. x 1-1/2" (3.81 cm)
- Syringe: 3 mL Luer-Lock
- Syringe: 10 ml with Locking NRFit™ Connector
- 10 mL Ampule 0.9% Saline Solution
- 5 mL Ampule 1% Lidocaine HCl Solution
- 5 mL Ampule 1.5% Lidocaine HCl with Epinephrine 1:200,000 Solution
- SharpsAway® Disposal Cup
- Towel: 13" x 18" (33 cm x 46 cm)
- Drape: 24" x 36" (61 cm x 91 cm) with 4" (10 cm) fenestration, with adhesive
- 5 Micron Particulate Filter Needle
- 5 Micron Particulate Filter Needle with NRFit™ connector
- Medication Label: 1% Lidocaine
- Medication Label: 0.9% Saline
- Medication Label: 1.5% Lidocaine with Epinephrine 1:200,000
- Gauze Pads: 4" x 4" (10 cm x 10 cm)



CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician. Contents of unopened, undamaged package are sterile. Disposable. Refer to package insert for current warnings, indications, contraindications, precautions, and Instructions for Use.

References:

1. Banwell B, et al. Decreased Incidence of Complications with the Arrow (FlexTip Plus) Epidural Catheter. *Canadian Journal of Anaesthesia*. 1998; 45: 370-72.
2. Cohen S, et al. Arrow versus B. Braun® Epidural Catheter for Epidural Block in Obstetric Patients. *Anesthesiology*. 2004; 101: A1199.
3. Cohen S, et al. Sutured Epidural Catheter for Epidural Block in Obstetric Patients: "Arrow" versus "B. Braun." *Anesthesiology*. 2006; 105: A904.
4. Juneja M, et al. Incidence of Epidural Vein Cannulation in Parturients with Three Different Epidural Catheters. *Regional Anesthesia & Pain Medicine*. 1996; 21 (2): 4.
5. Mhyre, Jill M, et al. A Systematic Review of Randomized Controlled Trials that Evaluate Strategies to Avoid Epidural Vein Cannulation during Obstetric Epidural Catheter Placement. *Anesthesia & Analgesia*. 2009; 108: 1232-1242.
6. Juneja M, et al. Comparison of Epidural Catheters Induced Paresthesia in Parturients. *Regional Anesthesia & Pain Medicine*. 1995; 20 (2): 152.

Teleflex is a global provider of medical technologies designed to improve the health and quality of people's lives. We apply purpose driven innovation – a relentless pursuit of identifying unmet clinical needs – to benefit patients and healthcare providers. Our portfolio is diverse, with solutions in the fields of vascular and interventional access, surgical, anesthesia, cardiac care, urology, emergency medicine and respiratory care. Teleflex employees worldwide are united in the understanding that what we do every day makes a difference. For more information, please visit teleflex.com.

Teleflex is the home of Arrow®, Deknatel®, Hudson RCI®, LMA®, Pilling®, Rüsçh®, and Weck® – trusted brands united by a common sense of purpose.

Corporate Office

Phone +1 610 225 6800, 550 E. Swedesford Road, Suite 400, Wayne, PA 19087, USA

Regional Offices

United States: Phone +1 919 544 8000, Toll Free 866 246 6990, cs@teleflex.com, 3015 Carrington Mill Boulevard, Morrisville, NC 27560, USA

Latin America: Phone +1 919 433 4999, la.cs@teleflex.com, 3015 Carrington Mill Boulevard, Morrisville, NC 27560, USA

International: Phone +353 (0)9 06 46 08 00, orders.intl@teleflex.com, Teleflex Medical Europe Ltd., IDA Business and Technology Park, Dublin Road, Athlone, Co Westmeath, Ireland

Australia +1300 360 226

Austria +43 (0)1 402 47 72

Belgium +32 (0)2 333 24 60

Canada +905 943 9004

China (Shanghai) +86 (0)21 6163 0965

China (Beijing) +86 (0)10 6418 5699

Czech Republic +420 (0)495 759 111

France +33 (0)5 62 18 79 40

Germany +49 (0)7151 406 0

Greece +30 210 67 77 717

India +91 (0)44 2836 5040

Italy +39 0362 58 911

Japan +81 (0)3 6632 3600

Korea +82 2 536 7550

Mexico +52 55 5002 3500

Netherlands +31 (0)88 00 215 00

New Zealand +0800 601 100

Portugal +351 22 541 90 85

Singapore (SEA non-direct sales countries) +65 6439 3000

Slovak Republic +421 (0)3377 254 28

South Africa +27 (0)11 807 4887

Spain +34 918 300 451

Switzerland +41 (0)31 818 40 90

United Kingdom +44 (0)1494 53 27 61

For more information, please visit teleflex.com.

NRFit™ is a trademark of GEDSA.

Teleflex, the Teleflex logo, Arrow, Deknatel, Hudson RCI, LMA, Pilling, Rüsçh, SnapLock, SharpsAway, TheraCath, and Weck are trademarks or registered trademarks of Teleflex Incorporated or its affiliates, in the U.S. and/or other countries. All other trademarks are trademarks of their respective owners.

Information in this document is not a substitute for the product Instructions for Use. The products in this document may not be available in all countries. Please contact your local representative. All data current at time of printing (07/2017). © 2017 Teleflex Incorporated. All rights reserved. MC-003322



Arrow®

Epidural Catheters with NRFit™ connectors

Enhancing patient comfort and safety

