

Angio Assist™

docking station

Angio Assist facilitates the introduction of guidewires into catheters and atherectomy burrs



teirstein EDGE™

device organizer

Teirstein Edge neatly organizes guidewires and catheters during interventional procedures




Organizing, securing and loading into

**Angio
Assist™**
docking station

User Friendly Design

- Friction-fit guidewire holder simplifies single-operator procedures
- Catheter loading area eliminates the need to touch or hold the stent during guidewire loading



Device Loading Areas

- Two slots facilitate alignment and introduction of 0.014" guidewires into catheters
- Three grooves facilitate alignment and introduction of 0.014" guidewires into 1.25mm–2.25mm atherectomy burrs

Unique S


- Designed with easy-to-use pi

Conventional tools has never been easier

teirstein
EDGE™
device organizer

Unique Grip Design

- Six friction fit slits secure guidewires and catheters
- Minimizes guidewire loss during procedure by tightly securing the wire in place
- Allows for controlled, micro movements of catheters and wires



Side Clips

Side clips for securing to drapes

Vessel Identification Labels

- Provide easy-to-identify guidewire and catheter location for individual vessels

Angio Assist™

docking station

The AngioAssist Docking Station is intended to stabilize and load guidewires and catheters during percutaneous procedures.

Ordering Specifications		
Model	Description	Box Quantity*
8151	Angio Assist docking station	16

*Individual units are not sold separately.



teirstein EDGE™

device organizer

The Teirstein Edge Device Organizer is intended to stabilize and load guidewires and catheters during percutaneous procedures.

Ordering Specifications		
Model	Description	Box Quantity*
8150	Teirstein Edge device organizer	10

*Individual units are not sold separately.



Please see the Instructions for Use for a complete listing of the indications, contraindications, warnings and precautions.

CAUTION: Federal Law (U.S.A.) restricts this device to sale by or on the order of a physician.

U.S. and International patents pending.

Angio Assist and Teirstein Edge are trademarks of Vascular Solutions, Inc.

©2013 Vascular Solutions, Inc. All rights reserved. ML2274 Rev. B 02/13



Vascular Solutions, Inc.
6464 Sycamore Court
Minneapolis, Minnesota 55369 USA

Customer Service:
United States: 888.240.6001
customerservice@vasc.com
www.vasc.com