

VSI Guidewires

0.018" and 0.014" Guidewires

Coil Tips
available in soft, flexible
stainless steel or highly
radiopaque tungsten

50cm 0.014"
for specialty use

Stainless Steel or Nitinol Mandrels
designed for flexibility and pushability

Separately Packaged
0.018" guidewires available
in a variety of lengths

VSI Guidewires

0.018" Guidewires

Stainless Steel Mandrel - Stainless Steel Soft Tip

Model	Wire Diameter	Wire Length	Special Features	Box Quantity
7300	0.018"	40cm	-	10
7301	0.018"	60cm	-	10

Nitinol Mandrel - Stainless Steel Soft Tip

Model	Wire Diameter	Wire Length	Special Features	Box Quantity
7302	0.018"	40cm	-	10
7303	0.018"	60cm	-	10

Stainless Steel Mandrel - Tungsten Soft Tip

Model	Wire Diameter	Wire Length	Special Features	Box Quantity
7304	0.018"	40cm	-	5
7305	0.018"	60cm	-	5
7326	0.018"	80cm	PTFE-coated	5
7327	0.018"	100cm	PTFE-coated	5

Nitinol Mandrel - Tungsten Soft Tip

Model	Wire Diameter	Wire Length	Special Features	Box Quantity
7306	0.018"	40cm	Angled Tip	5
7342	0.018"	40cm	-	5
7307	0.018"	60cm	-	5
7324	0.018"	80cm	-	5
7325	0.018"	130cm	Marked*	5

Devices packaged in the quantities listed above. Individual units are not sold separately.

*Reference markers at 5cm intervals between distal 20cm and 65cm

0.014" Guidewires

Nitinol Mandrel - Stainless Steel Soft Tip

Model	Wire Diameter	Wire Length	Special Features	Box Quantity
7343	0.014"	50cm	-	10

Packaged in quantities of 10 units per box. Individual units are not sold separately.

The VSI Guidewire is indicated for percutaneous entry of peripheral vessels using the Seldinger Technique. This device is not indicated for use in the coronary or cerebral vasculature. Please see the Instructions for Use for a complete listing of the indications, contraindications, warnings and precautions. CAUTION: Federal law (U.S.A.) restricts this device to sale by or on the order of a physician.

©2017 Vascular Solutions, Inc. All rights reserved. ML2440 Rev. D 02/17

This literature is intended for distribution only in the U.S.



Manufactured for:
Vascular Solutions, Inc.
 6464 Sycamore Court North
 Minneapolis, Minnesota 55369
 Customer Service: 888.240.6001
 www.vasc.com