


# OptiSeal<sup>®</sup>

## Valved Peelable Introducer

Optimal insertion and sealing during lead placement



Biphasic (two seal) valve designed to reduce risk of air embolism<sup>1</sup>

Symmetric, ergonomically designed wings allow for a smooth, unidirectional crack and peel

### Advanced safety & simplicity


- Designed to limit the clinician's exposure to blood-borne pathogens<sup>1</sup>
- Smooth sheath-to-dilator transition facilitates atraumatic vessel access
- Color coded by French size for easy identification:



<sup>1</sup>Vesely TM, Fazzaro AG, Gheradini D, Preliminary Evaluation of Valved Introducer Sheath for the Insertion of Tunneled Hemodialysis Catheters, *Seminars in Dialysis*, 2004; 17(1): 65-68

# OptiSeal® Valved Peelable Introducer

The OptiSeal Valved PTFE Peelable Introducer is intended for use in percutaneous insertion of pacing leads or catheters in the venous system.




ORDERING SPECIFICATIONS <sup>1</sup>					
Model	Size	Sheath I.D.	Sheath Length	Color	Box Quantity*
1003975-001	6F	0.086"	13cm	Dark Green	5
1003975-002	7F	0.099"	13cm	Orange	5
1003975-003	8F	0.112"	13cm	Light Blue	5
1003975-004	9F	0.125"	13cm	Black	5
1003975-005	10F	0.138"	13cm	White	5
1003975-006	11F	0.151"	13cm	Blue	5
1003975-007	12F	0.164"	13cm	Beige	5
1003975-008	6F	0.086"	25cm	Dark Green	5
1003975-009	7F	0.099"	25cm	Orange	5
1003975-010	8F	0.112"	25cm	Light Blue	5
1003975-011	9F	0.125"	25cm	Black	5
1003975-012	10F	0.138"	25cm	White	5
1003975-013	11F	0.151"	25cm	Blue	5
1003975-014	12F	0.164"	25cm	Beige	5

\*Individual units are not sold separately.

## Kit contents:

- 13cm or 25cm sheath & corresponding dilator
- 0.035", 60cm guidewire
- One 18G needle
- 12cc syringe



Please see the Instructions for Use for a complete listing of the indications, contraindications, warnings and precautions.

CAUTION: Federal Law (U.S.A.) restricts this device to sale by or on the order of a physician.

OptiSeal is a registered trademark of Greatbatch Ltd.

©2013 Vascular Solutions, Inc. All rights reserved. ML2255 Rev. D 08/13

This literature is intended for distribution only in the U.S.



**Distributed by Vascular Solutions, Inc.**  
6464 Sycamore Court North  
Minneapolis, Minnesota 55369 USA

Customer Service:  
United States: 888.240.6001  
customerservice@vasc.com  
www.vasc.com