



Hudson RCI®
ClearConnect™
 Secure Tubing Connector

Teleflex is improving respiratory care with every breath.

Inspired by physicians and clinicians who help develop products, we are proud to introduce Hudson RCI® ClearConnect™ Secure Tubing Connector, a new connector for oxygen therapy and aerosol therapy products.

- Threads onto the flowmeter directly and securely^{1,2}
- Soft interior allows it to be pushed onto barbed connections³
 - Pre-filled bubble humidifiers
 - Refillable bubble humidifiers
 - Barbed adaptors
 - Transport tanks
 - Compressors
- Clear adaptor nut in accordance with gas source guidelines^{4,5}
- Available on oxygen therapy and aerosol therapy products



Benefits



The Clinician
 Provides a new level of ease-of use and security to gas source connections for respiratory care^{1,2}



Your Institution
 Attached to a variety of oxygen and aerosol products and works on threaded and barbed connections



The Patient
 Reduces unwanted disconnections from flow meter gas source “pop-offs” and leaks less than a standard and universal connector^{1,2}

New connector technology for respiratory care



Standard Connector and Universal Connector compared to ClearConnect™ Secure Tubing Connector



1041CC



1115CC



1103CC



1883CC

Hudson RCI® ClearConnect™ Secure Tubing Connector

ITEM CODE	DESCRIPTION	PRODUCT LINE	CASE QTY
1041CC	Mask, medium concentration, adult elongated, 7' with CC	Oxygen Masks	50
1059CC	Mask, non-rebreather with safety vent, adult elongated, 7' with CC	Oxygen Masks	50
1115CC	Star Lumen® Tubing, 7' with CC	Oxygen Tubing	50
1103CC	Cannula, over-the-ear, 7' with CC	Nasal Cannulas	50
1810CC	Cannula, over-the-ear, 14' with CC	Nasal Cannulas	50
1883CC	Micro Mist® Neb Kit with mouthpiece, tee, reservoir tube, and 7' Star Lumen® Tubing	Small Volume Nebulizers	50
1734CC	Up-Draft II® Neb Kit with mouthpiece, tee, reservoir tube, and 7' Star Lumen® Tubing	Small Volume Nebulizers	50

References:

1. Data on file (2018 internal study), Teleflex Incorporated, Report #D023372 ClearConnect leak test comparison report. Bench testing may not be indicative of clinical performance.
2. Data on file (2018 internal study), Teleflex Incorporated, Report #D023373 ClearConnect flowmeter retention report. Bench testing may not be indicative of clinical performance.
3. Data on file (2018 internal study), Teleflex Incorporated, Report #D022731 ClearConnect push-on and ball flow meter validation testing report
4. Patient Safety and Usability: Experiences in the U.S. Department of Veterans Affairs. User Experience Magazine. http://uxpamagazine.org/patient_safety_usability/ Accessed March 26, 2018
5. Fatal gas line mix-up: How to avoid making this "gastly" mistake. Institute for Safe Medication Practices. http://www.ismp.org/Newsletters/acutecare/archives/ismpmsa12_16-04acutecare.pdf Accessed March 26, 2018
6. Patents: US 8,707,950; US D674,895; CA 2,806,056. Other patents pending.

Teleflex, the Teleflex logo, ClearConnect, Hudson RCI, Micro Mist, Star Lumen, and Up-Draft II are trademarks or registered trademarks of Teleflex Incorporated or its affiliates, in the U.S. and/or other countries. Information in this document is not a substitute for the product Instructions for Use. The products in this document may not be available in all countries. Please contact your local representative. All data current at time of printing (07/2018). Subject to technical changes without further notice. © 2018 Teleflex Incorporated. All rights reserved. MC-004395