



FLOCATH® AND FLOCATH QUICK™ HYDROPHILIC INTERMITTENT CATHETERS

Your solution for quick, easy intermittent catheterization

AT A GLANCE

- FloCath Quick Hydrophilic Intermittent Catheter with integrated 0.9% saline pouch and grip sleeve for easy activation, handling and touchless insertion
- Full range of sizes
- Available in Straight, Female and Coude
- Color-coded funnel to easily identify size
- Not made with natural rubber latex
- Compact, discrete packaging, complete with integrated saline pouch and thumb grip to ease opening

FLOCATH HYDROPHILIC INTERMITTENT CATHETER

The FloCath® Hydrophilic Intermittent Catheter features smooth, polished Soft-Eye technology and a gently tapered tip for comfortable insertion. The hydrophilic coating helps to reduce friction and tissue irritation that can occur with the use of uncoated catheters. The hydrophilic coating is easily activated by soaking the catheter in water for 30 seconds, eliminating possible contamination when applying lubricant to uncoated catheters.

FLOCATH QUICK HYDROPHILIC INTERMITTENT CATHETER

The RüsCh® FloCath Quick Catheter is an “All-In-One” intermittent catheter providing the ultimate convenience. The integrated package contains a FloCath Hydrophilic Catheter with our patented handling sheath and a packet of sterile 0.9% saline. Sterile saline is much

closer to the body’s natural fluids than sterile water, providing an extremely lubricious coating which allows for a much more comfortable insertion. The handling sheath provides an easy grip, allowing for touchless insertion. The sheath can also be pushed back and used as an extension to the catheter.

FLOCATH QUICK HYDROPHILIC INTERMITTENT CATHETER KITS

The FloCath Quick Kits offer the convenience of a small, discrete system while maintaining sterility at home or on the go. The compact kit includes a single, sterile FloCath Quick Catheter, Straight or Female, packaged complete with insertion supplies including gloves (not made with natural rubber latex), antiseptic skin prep, underpad, urine collection bag and a refuse bag.

ORDERING INFORMATION

FLOCATH® HYDROPHILIC INTERMITTENT CATHETER

FLOCATH HYDROPHILIC INTERMITTENT CATHETER, STRAIGHT TIP 30/BOX

PRODUCT #	SIZE
220800060	6 Fr.
220800080	8 Fr.
220800100	10 Fr.
220800120	12 Fr.
220800140	14 Fr.
220800160	16 Fr.
220900140	14 Fr. Female Length 7"

FLOCATH HYDROPHILIC INTERMITTENT CATHETER, COUDÉ TIP 30/BOX

PRODUCT #	SIZE
220850100	10 Fr.
220850120	12 Fr.
220850140	14 Fr.
220850160	16 Fr.

FLOCATH QUICK™ HYDROPHILIC INTERMITTENT CATHETER

FLOCATH QUICK HYDROPHILIC INTERMITTENT CATHETER, STRAIGHT 16" 30/BOX

PRODUCT #	SIZE
220400060	6 Fr.
220400080	8 Fr.
220400100	10 Fr.
220400120	12 Fr.
220400140	14 Fr.
220400160	16 Fr.

FLOCATH QUICK HYDROPHILIC INTERMITTENT CATHETER, FEMALE 7" 30/BOX

PRODUCT #	SIZE
220500060	6 Fr.
220500080	8 Fr.
220500100	10 Fr.
220500120	12 Fr.
220500140	14 Fr.
220500160	16 Fr.

FLOCATH QUICK HYDROPHILIC INTERMITTENT CATHETER, COUDÉ 16" 30/BOX

PRODUCT #	SIZE
220600100	10 Fr.
220600120	12 Fr.
220600140	14 Fr.
220600160	16 Fr.
220600180	18 Fr.

FLOCATH QUICK HYDROPHILIC INTERMITTENT CATHETER KIT

FLOCATH QUICK HYDROPHILIC INTERMITTENT CATHETER, STRAIGHT 16" 50/BOX

PRODUCT #	SIZE
221400060	6 Fr.
221400080	8 Fr.
221400100	10 Fr.
221400120	12 Fr.
221400140	14 Fr.
221400160	16 Fr.

FLOCATH QUICK HYDROPHILIC INTERMITTENT CATHETER, FEMALE 7" 50/BOX

PRODUCT #	SIZE
221500060	6 Fr.
221500080	8 Fr.
221500100	10 Fr.
221500120	12 Fr.
221500140	14 Fr.
221500160	16 Fr.

TO PLACE AN ORDER CONTACT CUSTOMER SERVICE AT 1.866.246.6990

Teleflex, FloCath, FloCath Quick and Rusch are trademarks or registered trademarks of Teleflex Incorporated or its affiliates. © 2014 Teleflex Incorporated. All rights reserved. 2014-2751

KIT CONTENTS

- Underpad
- One pair of large gloves not made with natural rubber latex
- Antiseptic skin prep
- Metered collection bag
- Refuse bag