



THE RUSCH BELLY BAG® URINARY COLLECTION DEVICE

Your Life, Your Health, Your Choice®

AT A GLANCE

- Materials contain no natural rubber latex
- Anti-Reflux valve prevents urine return from bag
- 52" belt can be trimmed to fit
- 1,000 mL capacity

The Rusch Belly Bag Urinary Collection Device is a solution for those wearing an indwelling catheter as they recover from surgery, ambulatory patients and those that lead a more active lifestyle and prefer a more discrete urine drainage bag. It helps to prevent accidental catheter extraction, leg bruising, skin break down and abrasions.

The Belly Bag Collection Device is a urine bag designed to be worn at the waist. It is designed for either males or females who have an indwelling Foley or suprapubic catheter. The bag fastens around the waist with a woven belt and a quick release buckle. It has a 1,000 mL capacity and is not manufactured from materials containing natural rubber latex. An anti-reflux valve behind the catheter port prevents reflux urine flow. The Belly Bag Collection Device is emptied by the drain tube with a twist valve and can be worn 24 hours a day. The bag is worn under the clothing where a soft backing provides a comfortable fit.

ORDERING INFORMATION

PRODUCT CODE	DESCRIPTION	BOX QTY.
B1000	Standard Belly Bag Collection Device with twist drain	10
B1000CT	Standard Belly Bag Collection Device with twist drain and 24" extension drain tube	10
B1000P	Standard Belly Bag Collection Device with twist drain and sample port	10

TO PLACE AN ORDER CONTACT CUSTOMER SERVICE AT 1.866.246.6990

Teleflex, Rusch and Your Life, Your Health, Your Choice are trademarks or registered trademarks of Teleflex Incorporated or its affiliates.
© 2013 Teleflex Incorporated. All rights reserved. 2012-1524